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**GB** **KKD® KONDISIL - Condensation Type Silicone Impression System**

**KKD® KONDISIL V-1 Putty - Silicone Tray Material**  
High-viscous Silicone Tray Material for preliminary impressions. Also for bite registrations and for laboratory procedures (as a space maintainer, as an insulating material etc.) very suitable. V-1 Putty material achieves in a few minutes an optimal shore hardness value and simultaneously develops exceptional elasticity. - Colour: white

**V-1 Putty Mixing rate:**  
1 measuring spoon (21 g) V-1 Putty to 2,5 cm (0,5 g) paste accelerator  
Mixing time: 1 minute  
Working time: 2 minutes  
Setting time in mouth: 3 ½ minutes  
Complete setting time: 5 ½ minutes

**Physical dates:**  
Shore hardness immediately: 70  
Shore hardness after 2 hours: 79  
Max. strain in compression: 1,3%  
Memory after deformation: 98,6%  
Linear dimensional change: 0,2%

Test conducted at room temperature of 23° C and relative humidity of 50 % (+/- 5%).  
EN ISO 4823:2000 Type 0

**Presentation:** REF:  
Pck.: 1,6 kg with hardener 18000S

**Individuals:**  
Jar: 1,6 kg without hardener 17998  
Pail: 10 kg without hardener 18001  
(please order approx. 7-10 units hardener to it)  
35 ml paste accelerator (red) 18003

as model material we recommend to use, according to final use, class III (e.g. Topdent Hartgips, KKD Cal) and class IV (e.g. KKD Satin Roc, KKD Veradur 03) gypsum

**KKD® KONDISIL V-3 light - Silicone Impression Material**  
This low-viscous Silicone Impression Material for refined impressions has an excellent fluidity and a precise detail reproduction. It can be used with an impression tray or with a syringe. When used in combination with the V-1 putty material an optimal compound is guaranteed. - Colour: blue

**V-3 low-viscous Mixing rate:**  
8 cm (6 g) V-3 base material to 5 cm (1 g) paste accelerator  
Mixing time: 45 to 60 sec.  
Working time: 2 minutes  
Setting time in mouth: 3 ½ minutes  
Complete setting time: 5 ½ minutes

**Physical dates:**  
Shore hardness immediately: 33  
Shore hardness after 2 hours: 39  
Max. strain in compression: 12,2%  
Memory after deformation: 99,2%  
Linear dimensional change: 1,1%

Test conducted at room temperature of 23° C and relative humidity of 50 % (+/- 5%).  
EN ISO 4823 Type 3

**Presentation:** REF:  
Pkg.: 200 ml with hardener and mixing pad 18008S

**Individuals:**  
200 ml V-3 base material 18009  
35 ml paste accelerator (red) 18003

as model material we recommend to use, according to final use, class III (e.g. Topdent Hartgips, KKD Cal) and class IV (e.g. KKD Satin Roc, KKD Veradur 03) gypsum

- smooth consistency
- easy to mix and to handle
- precise reproductions of details and borders
- great tensile strength
- elastic
- dimensionally stable
- mint-scented
- long shelf life

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