

ONECURE™



O-Star

WIDE-SPECTRUM CURING LIGHT
ULTIMATE CURING EXPERT

- All-Metal Body
- 1 Sec Curing
- Built-in Light Meter
- Uniform and Collimated Beam
- 7 Modes in Total

WOODPECKER Curing Light


16 years of professional technology
38 national patents
Global sales over 1.5 million units
Pioneering Curing Light expert in China

SMALL BUT POWERFUL MAKING A BIG DIFFERENCE

Full-metal body, ergonomically-streamlined design,
without tiredness from long-time holding
108g main unit,
light but strong, anti-shock and durable

7 MODES TO FULLY MEET A VARIETY OF CLINICAL NEEDS

Like an all star, O-Star is strong in performance, well-controlled and can fully meet various clinical needs.
Normal High Turbo Ortho Soft Pulse Check



Normal Mode
Constant light intensity: 1000-1200mw/cm²
Suitable for most treatment scenarios, such as fillings, restorations, etc.
Time setting: 5s, 10s, 15s, 20s

High Mode
Constant light intensity: 1800-2000 mw/cm²
Suitable for quick curing, ensuring enough curing depth.
Time setting: 3s, 5s

Turbo Mode
Time setting: 1s, 3s
2700-3000mw/cm²
Ultra-high light intensity enables 1 sec curing of 2mm resin.

Ortho Mode
Time setting: 3x5, 3x10
2700-3000mw/cm²
High light intensity with strong penetrability, especially suitable for orthodontic bracket bonding.

Soft Mode
Light intensity is gradually increased from 0 to 1200mw/cm²
Time setting: 5s, 10s, 15s, 20s
Effectively reduce the shrinkage rate of resin, Lower the risk of microleakage.

Pulse Mode
The LED will stop for 0.4s at an interval of 0.4s curing.
When working in a cycle of treatment, it can effectively reduce heat generation, efficiently dissipate heat and ensure the comfort of diagnosis and treatment.
Time setting: 5s, 10s, 15s, 20s

Check Mode
This mode outputs a single purple light. When special glasses are worn, the curies are obviously orange-red.
Application: Detection of caries, calculus and cracked tooth.

UNIFORM BEAM, STRONG PENETRABILITY


More-focused light output thanks to light angle of only 15.8°, ensuring good deep curing effect.

Curing distance (cm)	Light source	Light spot length (cm)	Light spot width (cm)	Light spot area (mm²)	Light spot area ratio	Light angle
5	Woodpecker	38	38	1444	-14.1%	15.80°
	VALO Cordless	41	41	1681		17.38°
	Power Cure	87	87	5941.665	253.5%	37.95°



LESS LIGHT LOSS

The light intensity of Curing Light will be reduced when the curing distance increases. Compared to other Curing Lights, O-Star undergoes less light loss.



CHARGING BASE WITH BUILT-IN LIGHT METER


There is a light meter built in the charging base. Accurate light intensity is required for efficient polymerization of the resin, bringing you and your patients the best treatment experience.



SIMPLE TWO-BUTTON DESIGN, ENABLES EASY OPERATION.


360° ROTATABLE HEAD MEETING VARIOUS CLINICAL NEEDS

360° rotatable head, easy to illuminate all angles in the mouth
The low-profile head can easily reach the posterior teeth.
The patient only needs to open the mouth slightly and enjoy comfortable treatment.




WIDE SPECTRUM CURING LIGHT CAN CURE ALL RESIN MATERIALS ON THE MARKET.

Camphorquinone (CQ) is used as photoinitiator in most of the resins on the market, but some resins use two kinds of new photoinitiators, i.e. TPO or Ivocerin. The main absorption wavelength of these three kinds of resins is 385-515nm. With a wavelength range of 385-515nm, O-Star is suitable for effective curing of most resins on the market.



UNPRECEDENTED CURING LIGHT INCORPORATING VARIOUS STRENGTHS



LIGHT INTENSITY UP TO 3000 MW/cm²

The high light intensity brings sufficient curing depth, which can reduce the curing time and greatly improve the efficiency of resin curing.

Equipped with 10wLED
Light intensity up to 3000 mw/cm², bringing you full confidence in curing

