Cercon ceram kiss Paste Liner

Cercon ceram love Paste Liner

**Ducera AllCeram Paste Liner** 

**Ducera press Paste Opaque** 

**Duceragold Kiss Paste Opaque** 

**Duceragold love Pastebonder Degunorm** 

**Duceragold love Paste Opaque** 

**Duceram Kiss Neutral Paste** 

**Duceram Kiss Paste Opaque** 

**Duceram love NE - Bonder** 

**Duceram love Paste Opaque** 

**Duceram Plus Neutral Paste** 

**Duceram Plus Paste Opaque** 

SYMBIOceram Paste Opaque

1.20 / GB Version: 23.02.2023 Revision date: Issue date: 31.10.2001 replaces version: 1.19

Page: 1/12

5360990231 Dentsply 130911 01954363

# SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. **Product identifier**

Trade name Cercon ceram kiss Paste Liner

> Cercon ceram love Paste Liner Ducera AllCeram Paste Liner Ducera press Paste Opaque Duceragold Kiss Paste Opaque

Duceragold love Pastebonder Degunorm

Material no.

Specification

VA-Nr

Duceragold love Paste Opaque **Duceram Kiss Neutral Paste** Duceram Kiss Paste Opaque Duceram love NE - Bonder Duceram love Paste Opaque **Duceram Plus Neutral Paste** Duceram Plus Paste Opaque SYMBIOceram Paste Opaque

REACH Registration No.: if available listed in Chapter. 3

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

Relevant applications identified For dental use only.

#### 1.3. Details of the supplier of the safety data sheet

Company DeguDent GmbH

> Postfach 1364 D-63403 Hanau

Telephone +49 (0)6181/59-5576 Telefax +49 (0)6181/59-5751

Email address SDB.Degudent-DE@dentsplysirona.com

#### 1.4. **Emergency telephone number**

**Emergency information** +49 (0)6181/59-50 (This telephone number is available during office

hours only.)

# SECTION 2: Hazards identification

#### 2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

**Cercon ceram kiss Paste Liner** 

Cercon ceram love Paste Liner

**Ducera AllCeram Paste Liner** 

**Ducera press Paste Opaque** 

**Duceragold Kiss Paste Opaque** 

**Duceragold love Pastebonder Degunorm** 

**Duceragold love Paste Opaque** 

**Duceram Kiss Neutral Paste** 

**Duceram Kiss Paste Opaque** 

**Duceram love NE - Bonder** 

Duceram Iove Paste Opaque

**Duceram Plus Neutral Paste** 

Ducerain Flus Neutral Paste

**Duceram Plus Paste Opaque SYMBIOceram Paste Opaque** 

 Version:
 1.20 / GB

 Revision date:
 23.02.2023

 Issue date:
 31.10.2001

 replaces version:
 1.19

 Page:
 2 / 12

 Material no.
 5360990231

 Specification
 130911

 VA-Nr
 01954363



Not a hazardous mixture according to Regulation (EC) No 1272/2008.

### 2.2. Label elements

# Labelling as per (EU) 1272/2008

Statutory basis Labelling not required according to EU-CLP Ordinance (1272/2008).

# 2.3. Other hazards

In case product dust is released:, Risk of infection of the lung after prolonged inhalation of dust particles.

A PBT/vPvB evaluation is not available, since a chemical safety evaluation is not required / has not been carried out.

# **SECTION 3: Composition/information on ingredients**

### **Chemical nature**

The preparation contains:, Vitreous porcelain, Pigments

### 3.1. Substances

-

### 3.2. Mixtures

Information on ingredients / Hazardous components as per EU-CLP Regulation (EC) No. 1272/2008

• butane-1,4-diol			< 15%
CAS-No.	110-63-4	EC-No.	203-786-5
Hexane-1,2,6-triol			> 15%
CAS-No.	106-69-4	EC-No.	203-424-6

Texts of H phrases, see in Chapter 16

### **SECTION 4: First aid measures**

# 4.1. Description of first aid measures

Take off all contaminated clothing immediately.

# Inhalation

Provide fresh air, seek medical advice if necessary.

### Skin contact

Wash off with plenty of water and soap, seek medical advice if necessary.

**Cercon ceram kiss Paste Liner** 

**Cercon ceram love Paste Liner** 

**Ducera AllCeram Paste Liner** 

**Ducera press Paste Opaque** 

**Duceragold Kiss Paste Opaque** 

**Duceragold love Pastebonder Degunorm** 

**Duceragold love Paste Opaque** 

**Duceram Kiss Neutral Paste** 

**Duceram Kiss Paste Opaque** 

**Duceram love NE - Bonder** 

**Duceram love Paste Opaque** 

Duceram Plus Neutral Paste

**Duceram Plus Paste Opaque** 

SYMBIOceram Paste Opaque

 Version:
 1.20 / GB

 Revision date:
 23.02.2023

 Issue date:
 31.10.2001

 replaces version:
 1.19

 Page:
 3 / 12

 Material no.
 5360990231

 Specification
 130911

 VA-Nr
 01954363



### Eye contact

With eye held open, thoroughly rinse immediately with plenty of water for at least 10 minutes. Seek medical advice.

### Ingestion

Clean mouth with water and drink afterwards plenty of water.

Call a physician immediately.

# 4.2. Most important symptoms and effects, both acute and delayed

### **Symptoms**

None known

### **Hazards**

None known

### 4.3. Indication of any immediate medical attention and special treatment needed

Symptomatic treatment.

# **SECTION 5: Firefighting measures**

# 5.1. Extinguishing media

Suitable extinguishing media: dry powder

CO<sub>2</sub>

Alcohol-resistant foam

water spray

# 5.2. Special hazards arising from the substance or mixture

Hazard-determining flue gases might develop in case of fire:

Carbon monoxide carbon dioxide

# 5.3. Advice for firefighters

Fire residues should be disposed of in accordance with the regulations.

In the event of fire, wear self-contained breathing apparatus.

# **SECTION 6: Accidental release measures**

### 6.1. Personal precautions, protective equipment and emergency procedures

Wear personal protective equipment.

Handle in accordance with good industrial hygiene and safety practice.

# 6.2. Environmental precautions

Do not allow the product into the following compartments:, sewage water, stretches of water, drainage systems

# 6.3. Methods and material for containment and cleaning up

Cercon ceram kiss Paste Liner

Cercon ceram love Paste Liner

**Ducera AllCeram Paste Liner** 

**Ducera press Paste Opaque** 

**Duceragold Kiss Paste Opaque** 

**Duceragold love Pastebonder Degunorm** 

**Duceragold love Paste Opaque** 

**Duceram Kiss Neutral Paste** 

**Duceram Kiss Paste Opaque** 

**Duceram love NE - Bonder** 

**Duceram love Paste Opaque** 

**Duceram Plus Neutral Paste** 

**Duceram Plus Paste Opaque** 

SYMBIOceram Paste Opaque

Version: Revision date: 23.02.2023 Issue date: 31.10.2001 replaces version: 1.19 Page: 4/12

Material no. 5360990231 Specification 130911 VA-Nr 01954363



Absorb with liquid-binding material, e.g. inert absorbent, sand, universal binding agents. Dispose of absorbed material in accordance with the regulations.

#### 6.4. Reference to other sections

Wear personal protective equipment; see section 8.

Disposal considerations; see section 13.

# **SECTION 7: Handling and storage**

#### 7.1. Precautions for safe handling

Avoid contact with skin and eyes.

Handle in accordance with good industrial hygiene and safety practice.

#### 7.2. Conditions for safe storage, including any incompatibilities

### Advice on protection against fire and explosion

Normal measures for preventive fire protection.

# Storage

Store in original container.

Keep tightly closed in a dry, cool and well-ventilated place.

Do not store together with: oxidizing agents

# German storage class

13 - Non Combustible Solids

#### 7.3. Specific end use(s)

We are unaware of any specific end uses which go beyond the data reported in Section 1.

# SECTION 8: Exposure controls/personal protection

#### 8.1. **Control parameters**

#### 8.2. **Exposure controls**

# **Engineering measures**

Ensure suitable suction/aeration at the work place and with operational machinery.

# Personal protective equipment

# Hand protection

Wear protective gloves made of the following materials:.

Polychloroprene (PCP) Glove material

Material thickness 0.5 mm Break through time > 480 min Glove material Nitrile rubb

Nitrile rubber/Nitrile latex (NBR)

Material thickness 0.35 mm Break through time > 480 min Glove material butyl-rubber

Cercon ceram kiss Paste Liner

Cercon ceram love Paste Liner

**Ducera AllCeram Paste Liner** 

**Ducera press Paste Opaque** 

**Duceragold Kiss Paste Opaque** 

**Duceragold love Pastebonder Degunorm** 

**Duceragold love Paste Opaque** 

**Duceram Kiss Neutral Paste** 

**Duceram Kiss Paste Opaque** 

**Duceram love NE - Bonder** 

**Duceram love Paste Opaque** 

**Duceram Plus Neutral Paste** 

**Duceram Plus Paste Opaque** 

SYMBIOceram Paste Opaque

1.20 / GB Version: 23.02.2023 Revision date: Issue date: 31.10.2001 replaces version: 1.19 Page: 5/12

Material no. 5360990231 Specification 130911 VA-Nr 01954363



Material thickness 0.5 mm Break through time > 480 min

Fluorinated rubber (FKM) Glove material

Material thickness 0.4 mm Break through time > 480 min

Polyvinyl chloride (PVC) Glove material

Material thickness 0.5 mm Break through time > 480 min

Source: GESTIS substance database (hazardous substance information system of commercial

professional associations)

# Eve/face protection

Safety glasses with side-shields

# Hygiene measures

Do not eat, drink, smoke, or sniff while at work. Wash your hands and/or face before breaks and before termination of work., Immediately change moistened and saturated work clothes.

# **Protective measures**

Handle in accordance with good industrial hygiene and safety practice.. If there is the possibility of skin/eye contact, the indicated hand/eye/body protection should be used.

# **SECTION 9: Physical and chemical properties**

#### Information on basic physical and chemical properties 9.1.

**Appearance** 

Form paste

Colour depending on staining result

Odour odourless

Odour threshold: not applicable

(40 g I)(20 °C) Ha ca. 10

Melting point/range not applicable

Preparation

Boiling point/range not applicable

Preparation

> 134 °C Flash point

Components with the lowest value

Information taken from reference works and the literature.

Evaporation rate not applicable, Preparation

**Cercon ceram kiss Paste Liner** 

**Cercon ceram love Paste Liner** 

Ducera AllCeram Paste Liner

**Ducera press Paste Opaque** 

**Duceragold Kiss Paste Opaque** 

**Duceragold love Pastebonder Degunorm** 

**Duceragold love Paste Opaque** 

**Duceram Kiss Neutral Paste** 

**Duceram Kiss Paste Opaque** 

**Duceram love NE - Bonder** 

**Duceram love Paste Opaque** 

**Duceram Plus Neutral Paste** 

**Duceram Plus Paste Opaque** 

SYMBIOceram Paste Opaque

 Version:
 1.20 / GB

 Revision date:
 23.02.2023

 Issue date:
 31.10.2001

 replaces version:
 1.19

 Page:
 6 / 12

 Material no.
 5360990231

 Specification
 130911

 VA-Nr
 01954363



Flammability (solid, gas) This product has not been classified as flammable.

Lower explosion limit not applicable

Upper explosion limit not applicable

Vapour pressure not applicable

Preparation

Density ca. 1.9 g/cm3 (20 °C)

Relative density no data available

Water solubility partly soluble

Partition coefficient: n-

octanol/water

no data available

Autoinflammability

Not capable of spontaneous combustion or heating.

Thermal decomposition no data available

Viscosity, dynamic not applicable

Preparation

Viscosity, kinematic no data available

Explosiveness no data available

Oxidizing properties no data available

9.2. Other information

Other information No further physicochemical data were determined.

# **SECTION 10: Stability and reactivity**

# 10.1. Reactivity

no data available

# 10.2. Chemical stability

Under normal conditions: stable.

### 10.3. Possibility of hazardous reactions

Possibility of hazardous Reaction with strong oxidants.

**Cercon ceram kiss Paste Liner** 

**Cercon ceram love Paste Liner** 

Ducera AllCeram Paste Liner

**Ducera press Paste Opaque** 

**Duceragold Kiss Paste Opaque** 

**Duceragold love Pastebonder Degunorm** 

**Duceragold love Paste Opaque** 

**Duceram Kiss Neutral Paste** 

**Duceram Kiss Paste Opaque** 

**Duceram love NE - Bonder** 

Duceram Iove Paste Opaque

**Duceram Plus Neutral Paste** 

**Duceram Plus Paste Opaque** 

**SYMBIOceram Paste Opaque** 

 Version:
 1.20 / GB

 Revision date:
 23.02.2023

 Issue date:
 31.10.2001

 replaces version:
 1.19

Page: 7 / 12

 Material no.
 5360990231

 Specification
 130911

 VA-Nr
 01954363



reactions

### 10.4. Conditions to avoid

None known.

# 10.5. Incompatible materials

None known

### 10.6. Hazardous decomposition products

No hazardous decomposition products known.

No further information available

# **SECTION 11: Toxicological information**

# 11.1. Information on toxicological effects

No results of animal experiments with the product available.

Acute oral toxicity LD50 rat: 1525 mg/kg

Test substance: butane-1,4-diol

Acute inhalation toxicity LC50 Rat: > 5.1 mg/l / 4 h

Test substance: butane-1,4-diol

Acute dermal toxicity LD50 Rabbit: > 2000 mg/kg

Test substance: butane-1,4-diol

Skin irritation Rabbit

not irritating

Method: OECD Test Guideline 404

Test substance: Hexane-1,2,6-triol

Eye irritation Rabbit

not irritating

Method: OECD Test Guideline 405

Test substance: Hexane-1,2,6-triol

Sensitization no data available

Repeated dose toxicity Oral Rat

Testing period: 90 d

NOEL: ca. 1000 mg/kg

target organ/effect: body weight development negative

Test substance: Hexane-1,2,6-triol

Feeding experiments

Oral Rat

Cercon ceram kiss Paste Liner

Cercon ceram love Paste Liner

Ducera AllCeram Paste Liner

**Ducera press Paste Opaque** 

**Duceragold Kiss Paste Opaque** 

**Duceragold love Pastebonder Degunorm** 

**Duceragold love Paste Opaque** 

**Duceram Kiss Neutral Paste** 

**Duceram Kiss Paste Opaque** 

**Duceram love NE - Bonder** 

**Duceram love Paste Opaque** 

Duceram Plus Neutral Paste

**Duceram Plus Paste Opaque** 

**SYMBIOceram Paste Opaque** 

 Version:
 1.20 / GB

 Revision date:
 23.02.2023

 Issue date:
 31.10.2001

 replaces version:
 1.19

Material no. 5360990231 Specification 130911 VA-Nr 01954363

**Dentsply** Sirona

Testing period: 2 Jahre

NOEL: 100 mg/kg
target organ/effect: Kidney, Liver

Test substance: Hexane-1,2,6-triol

Feeding experiments, chronic

Risk of aspiration toxicity 
Due to the viscosity, this product does not present an aspiration hazard.

Mutagenicity assessment no data available

8 / 12

Carcinogenicity No data available

Toxicity to reproduction No data available

Human experience Health injuries are not known or expected under normal use.

Further information No hazardous reactions are known if properly handled and stored., Risk of

infection of the lung after prolonged inhalation of dust particles.

# **SECTION 12: Ecological information**

### 12.1. Toxicity

Page:

No ecotoxicological data is available for this product.

Toxicity to fish LC50 Cyprinus carpio: > 1240.00 mg/l / 96 h

Test substance: butane-1,4-diol
Method: EC 92/69

LC0 Cyprinus carpio: 1240 mg/l / 96 h

Test substance: butane-1,4-diol Method: EC 92/69

Toxicity in aquatic

EC50 Daphnia magna: 813 mg/l / 48 h

invertebrates Test substance: butane-1,4-diol

Method: EC 92/69

NOEC Daphnia magna: >= 100.00 mg/l / 504 h

Test substance: butane-1,4-diol
Method: OECD 202 part 2

Toxicity to algae EC50 Desmodesmus subspicatus (green algae): > 1000.00 mg/l / 72 h

Test substance: butane-1,4-diol

**Cercon ceram kiss Paste Liner** 

**Cercon ceram love Paste Liner** 

Ducera AllCeram Paste Liner

**Ducera press Paste Opaque** 

**Duceragold Kiss Paste Opaque** 

**Duceragold love Pastebonder Degunorm** 

**Duceragold love Paste Opaque** 

**Duceram Kiss Neutral Paste** 

**Duceram Kiss Paste Opaque** 

**Duceram love NE - Bonder** 

Duceram love Paste Opaque

Duceram Plus Neutral Paste

Duceram Plus Paste Opaque

SYMBIOceram Paste Opaque

 Version:
 1.20 / GB

 Revision date:
 23.02.2023

 Issue date:
 31.10.2001

 replaces version:
 1.19

 Page:
 9 / 12

 Material no.
 5360990231

 Specification
 130911

 VA-Nr
 01954363

**Dentsply** Sirona

Method: EC 92/69

NOEC Desmodesmus subspicatus (green algae): >= 1000.00 mg/l / 72 h

Test substance: butane-1,4-diol
Method: EC 92/69

Toxicity to bacteria EC 10 Pseudomonas putida / 16 h

Test substance: butane-1,4-diol

Method: Bringmann und Kühn, Z. Wasser Abwasser Forsch. 10, 87-

98 (1977)

# 12.2. Persistence and degradability

Biodegradability Result: 96 %

Test substance: butane-1,4-diol
Method: OECD 301 E

readily biodegradable

### 12.3. Bioaccumulative potential

Bioaccumulation no data available

12.4. Mobility in soil

Mobility No data available

### 12.5. Results of PBT and vPvB assessment

A PBT/vPvB evaluation is not available, since a chemical safety evaluation is not required / has not been carried out.

### 12.6. Other adverse effects

Further Information No data can be given because of insolubility in water.

Do not allow entrance in sewage water and soil.

Cercon ceram kiss Paste Liner

Cercon ceram love Paste Liner

**Ducera AllCeram Paste Liner** 

**Ducera press Paste Opaque** 

**Duceragold Kiss Paste Opaque** 

**Duceragold love Pastebonder Degunorm** 

**Duceragold love Paste Opaque** 

**Duceram Kiss Neutral Paste** 

**Duceram Kiss Paste Opaque** 

**Duceram love NE - Bonder** 

**Duceram love Paste Opaque** 

**Duceram Plus Neutral Paste** 

Duceralli Flus Neutral Faste

**Duceram Plus Paste Opaque SYMBIOceram Paste Opaque** 

 Version:
 1.20 / GB

 Revision date:
 23.02.2023

 Issue date:
 31.10.2001

 replaces version:
 1.19

 Page:
 10 / 12

 Material no.
 5360990231

 Specification
 130911

 VA-Nr
 01954363



# **SECTION 13: Disposal considerations**

### 13.1. Waste treatment methods

### **Product**

Disposal according to local authority regulations.

# **Uncleaned packaging**

Disposal according to local authority regulations.

# **SECTION 14: Transport information**

Not dangerous according to transport regulations.

- 14.1. UN number: --
- 14.2. UN proper shipping name: --
- 14.3. Transport hazard class(es):
- 14.4. Packing group:
- 14.5. Environmental hazards: --
- 14.6 Special precautions for user: No

### **SECTION 15: Regulatory information**

# 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture National legislation

### 15.2. Chemical safety assessment

Chemical safety assessment No Chemical Safety Report as per Articles 2(8), 2(9) or 14 of the REACH Regulatione is required for this product.

# **SECTION 16: Other information**

Classification and applied procedure to derive the classification of mixtures according to EU Regulation (EC) No. 1272/2008 (CLP)

### **Further information**

Changes since the last version are highlighted in the margin. This version replaces all previous versions.

Cercon ceram kiss Paste Liner Cercon ceram love Paste Liner Ducera AllCeram Paste Liner Ducera press Paste Opaque

**Duceragold Kiss Paste Opaque** 

**Duceragold love Pastebonder Degunorm** 

Duceragold love Paste Opaque Duceram Kiss Neutral Paste Duceram Kiss Paste Opaque Duceram love NE - Bonder Duceram love Paste Opaque Duceram Plus Neutral Paste Duceram Plus Paste Opaque SYMBIOceram Paste Opaque

 Version:
 1.20 / GB

 Revision date:
 23.02.2023

 Issue date:
 31.10.2001

 replaces version:
 1.19

 Page:
 11 / 12

 Material no.
 5360990231

 Specification
 130911

 VA-Nr
 01954363



This information and all further technical advice is based on our present knowledge and experience. However, it implies no liability or other legal responsibility on our part, including with regard to existing third party intellectual property rights, especially patent rights. In particular, no warranty, whether express or implied, or guarantee of product properties in the legal sense is intended or implied. We reserve the right to make any changes according to technological progress or further developments. The customer is not released from the obligation to conduct careful inspection and testing of incoming goods. Performance of the product described herein should be verified by testing, which should be carried out only by qualified experts in the sole responsibility of a customer. Reference to trade names used by other companies is neither a recommendation, nor does it imply that similar products could not be used.

Legend

ADR European Agreement concerning the International Carriage of Dangerous Goods by

Road

ADN European Agreement concerning the International Carriage of Dangerous Goods by

**Inland Waterways** 

**ASTM** American Society for Testing and Materials

ATP Adaptation to Technical Progress

**BCF** Bioconcentration factor

BetrSichV German Ordinance on Industrial Safety and Health

**c.c.** closed cup

**CAS** Chemical Abstract Services

**CESIO** European Committee of Organic Surfactants and their Intermediates

ChemG German Chemicals Act

**CMR** carcinogenic-mutagenic-toxic for reproduction

**DIN** German Institute for Standardization

DMEL Derived minimum effect level
DNEL Derived no effect level

**EINECS** European Inventory of Existing Commercial Chemical Substances

**EC50** half maximal effective concentration

**GefStoffV** German Ordinance on Hazardous Substances

GGVSEB German ordinance for road, rail and inland waterway transportation of dangerous

goods

GGVSee German ordinance for sea transportation of dangerous goods

GLP Good Laboratory Practice
GMO Genetic Modified Organism

IATA International Air Transport Association
 ICAO International Civil Aviation Organization
 IMDG International Maritime Dangerous Goods
 ISO International Organization For Standardization

**LOAEL** Lowest observed adverse effect level

LOELLowest observed effect levelNOAELNo observed adverse effect levelNOECno observed effect concentration

**NOEL** no observed effect level

Cercon ceram kiss Paste Liner

**Cercon ceram love Paste Liner** 

**Ducera AllCeram Paste Liner** 

**Ducera press Paste Opaque** 

**Duceragold Kiss Paste Opaque** 

**Duceragold love Pastebonder Degunorm** 

**Duceragold love Paste Opaque** 

**Duceram Kiss Neutral Paste** 

**Duceram Kiss Paste Opaque** 

**Duceram love NE - Bonder** 

Duceram love Paste Opaque

Duceram Plus Neutral Paste

**Duceram Plus Paste Opaque** 

SYMBIOceram Paste Opaque

 Version:
 1.20 / GB

 Revision date:
 23.02.2023

 Issue date:
 31.10.2001

 replaces version:
 1.19

 Page:
 12 / 12

Material no. 5360990231 Specification 130911 VA-Nr 01954363



o. c. open cup

**OECD** Organisation for Economic Cooperation and Development

OEL Occupational Exposure Limit
PBT Persistent, bioaccumulative, toxic
PEC Predicted effect concentration
PNEC Predicted no effect concentration

**REACH** REACH registration

RID Convention concerning International Carriage by Rail

STOT Specific Target Organ Toxicity
SVHC Substances of Very High Concern

TA Technical Instructions

**TPR** Third Party Representative (Art. 4)

TRGS Technical Rules for Hazardous Substances
VCI German chemical industry association
vPvB very persistent, very bioaccumulative

**VOC** volatile organic compounds

VwVwS German Administrative Regulation on the Classification of Substances Hazardous to

Waters into Water Hazard Classes

WGK Water Hazard Class
WHO World Health Organization