

The perfect preventive dental care combination: KaVo SMARTmatic and Kerr attachments.



The new KaVo SMARTmatic
Prevention range.

- **Made in Germany** high product quality
- **Small head** for optimal visibility
- **Oscillating movement** for effective polishing
- **High-quality materials** for long life
- **Smart introduction** into the KaVo preventive healthcare world

* Applies only for PROPHY S33 and S53

Kerr cleaning and polishing accessories for preventive dental care.

- **Wide range** for maximum flexibility with preventive treatments
- **Innovative design** for the continuous application of polish onto teeth
- **Easy to use** for a smooth preventive process
- **Ideal** for the latest KaVo SMARTmatic dental instruments

KAVO  **Kerr**

The smart combination



SMARTmatic model	Prophy S19	Prophy S31	Prophy S33	S53
Speed transmission	8:1	8:1	8:1 / 70° osc. Rotary motion	
Max. Speed drive rpm.	20,000	20,000	20,000	20,000
Shaft/clamping mechanism	Doriot	Screw-In, Snap-On	Snap-On	Latch/2.35 mm
Item number	1.011.6740	1.011.6760	1.011.6800	1.011.6810

Selection of compatible Kerr attachments:

Pro-Cup® and Pro-Cup® Junior

✓ Compatible
— Not compatible

Latch type	Soft	Hard	Junior - Soft	Junior - Hard	Prophy S19	Prophy S31	Prophy S33	S53
					—	—	—	✓
Screw-type	Soft	Hard	Junior - Soft	Junior - Hard	Prophy S19	Prophy S31	Prophy S33	S53
					—	✓	—	—

Small dental and polishing brushes

Miniature Brushes	Natural	Nylon	Prophy S19	Prophy S31	Prophy S33	S53
			—	—	—	✓
Pro Brush	Natural	Nylon	Prophy S19	Prophy S31	Prophy S33	S53
			—	—	—	✓

Prophy cups

Snap-On	Ribs	Ribbed	Prophy S19	Prophy S31	Prophy S33	S53
			—	✓	✓	—

Prophy brushes

Snap-On	Natural bristles, normal	Natural bristles, pointed	Prophy S19	Prophy S31	Prophy S33	S53
			—	✓	✓	—



Scan and discover the complete Kerr portfolio: Attachments, cups and brushes for an efficient and splash-free polish application. (Costs in accordance with mobile phone agreement)

The products and equipment shown and described in this catalogue, as well as the services, are not available in every country. All information is accurate at the time of going to print. KaVo Dental GmbH assumes no responsibility whatsoever for variances in colour or shape from the depictions, nor for mistakes such as printing errors and reserves the right to make changes to the catalogue at any time. Reprinting, even of sections, is only allowed with permission from KaVo Dental GmbH.