# SECTION 1. Product identification/preparation/company undertaking

1.1 Product identifier:

Product no.: 81260

Trade name: BISICO Softbase Primer

Trade number: 13860

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against:

No further relevant information available.

#### Application of the substance / the mixture:

Dental-soft-reliner

### 1.3. Details of the supplier of the safety data sheet:

BISICO Bielefelder Dentalsilicone GmbH & Co. KG

Johanneswerkstraße 3 D-33611 Bielefeld Tel.: +49 521 8016800

Fax: +49 521 8016801 Email: info@Bisico.de

responsible for the Safety Data Sheet: Sebastian Zimmermann

#### 1.4. Emergency telephone number:

BISICO Bielefelder Dentalsilicone GmbH & Co. KG

Tel.: +49 521 8016800 (08:00-16:00 PM)

Fax: +49 521 8016801 S. Zimmermann

# **SECTION 2. <u>Hazards identification</u>**

### 2.1 Classification of the substance or mixture:

#### 2.1.1 Classification according to Regulation (EC) No 1272/2008:

Flam. Liq. 2 H225 Highly flammable liquid and vapour.

Skin Irrit. 2 H319 Causes serious eye irritation

Skin Sens. 1 H317 May cause an allergic skin reaction.

STOT SE 3 H336 May cause drowsiness or dizziness

# 2.2 Label elements:

**Labelling according to Regulation (EC) No 1272/2008:** The product is classified and labelled according to the CLP regulation.

#### Hazard pictograms:



GHS02



Signal word:

Danger

# bisico

# Safety Data Sheet according to 1907/2006/EC, Article 31

#### Hazard statements:

H225 Highly flammable liquid and vapour.

H319 Causes serious eye irritation

H336 May cause drowsiness or dizziness

### **Precautionary statements:**

P210 Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

#### Labeling according to EU guidelines:

Medical devices as defined in Directive 93/42/EEC and which are invasive or used in direct physical contact with the human body; are exempted from the provisions of Regulation (EC) No. 1272/2008 (CLP/GHS) usually if they are in the finished state and intended for the final user.

#### 2.3 Other hazards:

Results of PBT and vPvB assessment:

**PBT:** not applicable **vPvB:** not applicable

### SECTION 3. Composition/Information of Ingredients

#### 3.1 Substances:

3.2 Mixtures:

#### Chemical characterization:

#### **Description:**

Solution of polyacrylate in ethylacetate

#### **Hazardous components:**

CAS-No. [141-78-6],EC-No. 2005-500-4

Ethylacetate GHS02, GHS07, H225, H319, H336

#### **SECTION 4. First-aid measures**

# 4.1 Description of first aid measures:

general information: no special measures required

after inhalation: Supply fresh air

**after skin contact:** wash with water and

**after eye contact:** Rinse opened eye for several minutes under running water

and consult a docton

after swallowing: Caution if victim vomits. Risk of aspiration. Keep airways

free. Consult a physician. Subsequently administer: activated charcoal (20 – 40 g in 10% slurry). Laxative:

Sodium sulfate (1 tablespoon/ 1/4 I water).

### 4.2 Most important symptoms and effects, both acute and delayed:

No further relevant information available.

# 4.3 Indication of any immediate medical attention and special treatment needed:

No further relevant information available.

# **SECTION 5. Fire fighting measures**

#### 5.1 Extinguishing media:

**Suitable extinguishing agents:** CO<sub>2</sub>, powder, water spray. Use extinguishing method and media suitable to surrounding fire

**5.2 Special hazards arising from the substance or mixture:** No further relevant information available.

#### 5.3 Advice for firefighters:

# **Protective equipment:**

Wear self-contained respiratory protective device.

#### **SECTION 6. Accidental release measures**

### 6.1 Personal precautions, protective equipment and emergency procedures:

Do not inhale vapours/aerosols. Avoid substance contact. Ensure supply of fresh air in enclosed rooms

# 6.2 Environmental precautions:

Do not allow to enter drains/surface water/ ground-water.

#### 6.3 Methods for cleaning up:

Absorb with liquid binding material. Clean up affected area.

#### **6.4 Reference to other sections:** See Section 7 for information on safe handling

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

# **SECTION 7. Handling and storage**

# Handling:

This product should only be supplied or handled to dentists or dental technicians or upon their instructions.

#### 7.1 Precautions for safe handling:

Observe normal care for working with chemicals.

### Information about fire - and explosion protection:

No special measures required

#### 7.2 Conditions for safe storage, including any incompatibilities:

#### Requirements to be met by storerooms and receptacles:

Store only in the original containers

#### Further information about storage conditions:

Store in well ventilated place at 15 °C to 25 °C, protected from light. Store away from food stuffs. Store tightly closed. Keep away from heat and sources of ignition.

**7.3 Specific end use(s):**No further relevant information available.

# SECTION 8. Exposure controls/personal protection

#### 8.1. Control parameters:

Ingredients with limit values that require monitoring at the workplace:

CAS-No.	name	art		value ui	<u>nit</u>
[141-78-6]	Ethylacetate	TRGS900	400	ml/m <sup>3</sup>	
		or		1500 mg	$_{ m /m}^{ m 3}$

**Additional information:** The lists valid during the making were used as basis.

#### 8.2 Exposure controls:

#### Personal protective equipment

### General protective and hygienic measures:

The usual precautionary measures are to be adhered to when handling chemicals. Keep away from foodstuffs and beverages. Avoid contact with eyes or skin

# Hand protection:

Protective gloves.

Due to missing tests no recommendation to the glove material can be given for the product.

Eye protection: safety glasses

Respiratory protection: Required, when aerosols are generated

#### **SECTION 9. Physical and chemical properties**

#### 9.1. Information on basic physical and chemical properties:

**Appearance** 

Form: liquid
Color: slightly yellow
Odor: sweetish

Value Unit Method

Safety relevant basic data

Melting point/melting range: n.a.

**Boiling point/boiling range:** approx. 77 °C

Flash point: -4 °C DIN 51755

Ignition temperature: No information available

Danger of explosion:

Lower limit of explosion 2.1 Vol% Kühn-Upper limit of explosion 11.5 Vol% Birett ca. 0,9 (20°C) g/cm<sup>3</sup>

Vapor pressure: n.d.

Solubility in/miscibility with:

Water: soluble

**Content of solvents:** 

Organic solvents: ethylacetate

9.2 Other information No further relevant information available.

# SECTION 10. Stability and reactivity

10.1. Reactivity:

10.2. Chemical stability:

10.3 Possibility of hazardous reactions: no information available

#### 10.4 Conditions to avoid:

Heat, direct sunlight and radical starters. The product is stable by proper use. Avoid contact with alkali metals, hydrids, strong oxidizing agents, radical intitiators.

10.5 Incompatible materials:

**10.6 Hazardous decomposition products:** No information available.

### **SECTION 11.** Toxicological information

#### 11.1. Information on toxicological effects:

Acute toxicity:

 $LD_{50}$  (oral, rat) 5600 mg/kg  $LD_{50}$  (dermal, rabbit) >18000 mg/kg  $LD_{50}$  (inhal., rat) 5,86 mg/l (8h)

Primary irritant effect:

**Skin corrosion/irritation:** repeated exposure may cause skin dryness or cracking

Serious eye damage/irritation: strong irritating effect

**Inhalation and swallowing:** irritation of mucous membranes, loss of appetite, headache,

drowsiness

### Additional toxicological information:

Allergic reactions to acrylates after long-term exposure are possible. When used and handled according to specifications, the product does not have any harmful effects to our experience and the information provided to us.

#### **SECTION 12.** <u>Ecological information</u>

12.1 Toxicity:

**Acquatic toxicity:** No further relevant information available.

**12.2 Persistence and degradability:**No further relevant information available.

Behaviour in environmental systems:

**12.3 Bioaccumulative potential:** No further relevant information available.

**12.4 Mobility in soil:**No further relevant information available.

Additional ecological information:

General notes: Do not flush into surface water or sanitary sewer system

**12.5 Results of PBT and vPvB assessment: PBT:** Not applicable.

**vPvB:** Not applicable.

**12.6 Other adverse effects:** No further relevant information available.

#### **SECTION 13. Disposal consideration**

13.1. Waste treatment methods:

**Product:** 

Recommendation:

Use proper landfill disposal or incineration in accordance with local, state and federal

regulations.

**Uncleaned packaging:** 

**Recommendation:** Disposal must be made according to official regulations

#### **SECTION 14. Transport information**

14.1 UN-Number

ADR, IMDG, IATA: UN1173

14.2 UN proper shipping name

ADR:

**ETHYLACETAT** 

IMDG:

**ETHYLACETAT** 

IATA:

**ETHYLACETATE** 

14.3 Transport hazard class(es)

ADR:

UN 1173 ETHYLACETAT, 3, II

IMDG, IATA:

UN 1173 ETHYLACETAT, 3, II

Air transport ICAO-TI and IATA-DGR:

UN 1173 ETHYLACETATE, 3, II

14.4 Packing group ADR, IMDG, IATA:

Ш

14.5 Environmental hazards:

Marine pollutant: No

bisico

#### 14.6 Special precautions for user:

# 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code:

The transport regulations are cited according to international regulations and in the form applicable in Germany. Possible national deviations in other countries are not considered.

# **SECTION 15. Regulatory information**

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture:

#### Classification according to EC-guidelines:

The product is a medical device according to EC-directive 93/42 EEC

**15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

### **SECTION 16. Other information:**

Changes compared with the previous version: Adaptation according to 1907/2006/EG, Article 31

The above information is based on our present day knowledge and relates solely to the safety requirements of the product. The data do not signify any warranty with regards to products properties. However users of the product should satisfy themselves that the information given is sufficient and correct for their specific circumstances of use.

# Preparation by:

S. Zimmermann

Tel.: +49 521 8016800

1. <u>Product identification/preparation/company undertaking</u>

**Product details** 

Product no.: 81240

Trade name: BISICO Softbase Glanzlack

Trade number: 13850

Use or preparation: Dental-Soft-Reliner

Manufacturer/supplier: BISICO Bielefelder Dentalsilicone GmbH & Co. KG

Johanneswerkstraße 3 D-33611 Bielefeld

Emergency information: Tel.: +49 521 8016800

Fax: +49 521 8016801

2. <u>Hazards identification</u>

Classification of the substance or mixture:

Classification according to Regulation (EC)

1272/2008:

void

Information concerning particular hazards for human and environment:

The product does not have to be labeled due to the calculation procedure of the "General classification guideline for preparations of the EU" in the latest valid

version.

Label elements:

Labeling according to Regulation (EC)

1272/2008:

void

Labeling according to EU guidelines:

Medical devices as defined in Directive 93/42/EEC and which are invasive or used in direct physical contact with the human body; are exempted from the provisions of Regulation (EC) No. 1272/2008 (CLP/GHS) usually if they are in the finished state and intended for the final

user.

Other hazards: PBT: not applicable

Results of PBT and vPvB assessment: vPvB: not applicable

3. Composition/Information of Ingredients

**Chemical characterization:** 

**Description:** mixture of polydimethyl polymethyl vinyl siloxane,

polydimethyl polymethyl hydrogen siloxane and silica

Hazardous components: none

# bisico

4. First-aid measures

**General information:** no special measures required **after skin contact:** wash with plenty of water and soap

after eye contact: rinse with plenty of water and consult a doctor

after swallowing: If symptoms persist consult a doctor

5. Fire fighting measures

**Extinguishing media:** CO<sub>2</sub>, water

**Protective equipment:** no special measures required

6. <u>Accidental release measures</u>

Personal precautions: not required

**Environmental precautions:** no special measures required

Methods for cleaning up: pick up mechanically

Additional information: none

7. <u>Handling and storage</u>

Handling:

**Information for safe handling:** For dental use only. No special measures required

Recommendation for fire and

explosion protection:

No special measures required

Storage:

Requirements at storerooms and

containers:

No special measures required

Requirements for storage with other

products:

not required

Further information on storage

conditions:

none

8. <u>Exposure controls/personal protection</u>

Components with limits of values to be supervised at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at

the workplace

**Additional information:** The lists valid during the making were used as basis.

Personal protective equipment:

General measures of protection and

hygiene:

normal hygienic measures

not required Respiratory protection: not required Protection of hands: not absolute required Eye protection: 9. **Physical and chemical properties** Appearance: Paste Form: clear Color: odorless Odor: Value Unit Method Information on change in the physical °C Melting point/melting range: n.a. °C Boiling point/boiling range: n.a. °C n.a. Flash point: **Autoignition temperature:** n.a. Danger of explosion: none **Density:** ca. 1,0 (20°C) g/cm3 mbar Vapor pressure: n.a. Paste Viscosity: neutral pH: Partially soluble in toluene, petrol ether Solubility in/miscibility with: Water: insoluble **Content of solvents:** none **Organic solvents:** Water: Content of solids: n.a. 10. **Stability and reactivity** No decomposition if used according to specification Conditions to avoid: **Hazardous decomposition products:** none under normal conditions of storage and use 11. **Toxicological information** Acute toxicity: **Primary irritation:** 

no irritating effect skin: no irritating effect eye: Additional toxicological information: When used and handled according to specifications, the product does not have any harmful effects to our experience and the information provided to us. 12. **Ecological information General information:** Avoid transfer into environment WGK 1 (German regulation) slightly hazardous for water Classification of water endangerment: 13. **Disposal consideration Product:** Small quantities can be disposed of with household Recommendation: waste. Use proper landfill disposal or incineration in accordance with local, state and federal regulations **Uncleaned packaging:** Disposal must be made according to official regulations. Recommendation: 14. **Transport information** Land transport ADR, RID: Maritime transport: IMDG-Code: Air transport ICAO-TI/IATA-DGR: 15. Regulatory information The product is a medical device according to EC-Classification according to ECdirective 93/42 EEC. guidelines: A Chemical Safety Assessment has not been carried out. **Chemical safety assessment:** 16. **Other information:** Changes compared with the previous Adaptation according to 1907/2006/EG, Article 31 version: The above information is based on our present day knowledge and relates solely to the safety requirements of the product. The data do not signify any warranty with regards to products properties. However users of the product should satisfy themselves that the information given is sufficient and correct for their specific circumstances of use.

S. Zimmermann Tel.: +49 521 8016800

Preparation by:

1. <u>Product identification/preparation/company undertaking</u>

**Product details** 

Product no.: 81211

Trade name: BISICO Softbase cartridge

Trade number: 93810

Use or preparation: Dental-Soft-Reliner

Manufacturer/supplier: BISICO Bielefelder Dentalsilicone GmbH & Co. KG

Johanneswerkstraße 3 D-33611 Bielefeld

Emergency information: Tel.: +49 521 8016800

Fax: +49 521 8016801

2. <u>Hazards identification</u>

Classification of the substance or mixture:

Classification according to Regulation (EC)

1272/2008:

void

Information concerning particular hazards for human and environment:

The product does not have to be labeled due to the calculation procedure of the "General classification guideline for preparations of the EU" in the latest valid

version.

Label elements: -

Labeling according to Regulation (EC)

1272/2008:

**Labeling according to EU guidelines:** Medical devices as defined in Directive 93/42/EEC and

which are invasive or used in direct physical contact with the human body; are exempted from the provisions of Regulation (EC) No. 1272/2008 (CLP/GHS) usually if they are in the finished state and intended for the final

user.

Other hazards: PBT: not applicable

Results of PBT and vPvB assessment: vPvB: not applicable

3. Composition/Information of Ingredients

**Chemical characterization:** 

**Description:** mixture of polydimethyl polymethyl vinyl siloxane,

polydimethyl polymethyl hydrogen siloxane and silica

Hazardous components: none

bisico

4. First-aid measures

> General information: no special measures required

after skin contact: wash with plenty of water and soap

after eye contact: rinse with plenty of water and consult a doctor

after swallowing: If symptoms persist consult a doctor

5. Fire fighting measures

> **Extinguishing media:** CO<sub>2</sub>, water

**Protective equipment:** no special measures required

6. Accidental release measures

> Personal precautions: not required

**Environmental precautions:** no special measures required

pick up mechanically Methods for cleaning up:

Additional information: none

7. Handling and storage

Handling:

For dental use only. No special measures required Information for safe handling:

Recommendation for fire and

explosion protection:

No special measures required

Storage:

Requirements at storerooms and

containers:

No special measures required

Requirements for storage with other

products:

not required

Further information on storage

conditions:

none

8. **Exposure controls/personal protection** 

Components with limits of values to be

supervised at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at

the workplace

The lists valid during the making were used as basis. Additional information:

Personal protective equipment:

bisico

General measures of protection and

hygiene:

normal hygienic measures

Respiratory protection: not required

Protection of hands: not required

Eye protection: not absolute required

9. Physical and chemical properties

Appearance:

Form: Paste

Color: A-Comp.: white, B-Comp.: gum-colored

Odor: odorless

Value Unit Method

Information on change in the physical

state

Melting point/melting range: n.a. °C

Boiling point/boiling range: n.a. °C

Flash point: n.a. °C

**Autoignition temperature:** n.a.

Danger of explosion: none

**Density:** 1,05 (20°C) g/cm<sup>3</sup>

Vapor pressure: n.a. mbar

Viscosity: Paste

pH: neutral

Solubility in/miscibility with: Partially soluble in toluene, petrol ether

Water: insoluble

Content of solvents: none

Organic solvents: -

Water: -

Content of solids: n.a.

10. Stability and reactivity

Conditions to avoid: No decomposition if used according to specification

Hazardous decomposition products: none under normal conditions of storage and use

# bisico

11. <u>Toxicological information</u>

Acute toxicity: -

Primary irritation: -

skin: no irritating effect

eye: no irritating effect

Additional toxicological information: When used and handled according to specifications, the

product does not have any harmful effects to our experience and the information provided to us.

12. <u>Ecological information</u>

General information: Avoid transfer into environment

Classification of water endangerment: WGK 1 (German regulation) slightly hazardous for water

13. <u>Disposal consideration</u>

**Product:** 

**Recommendation:** Small quantities can be disposed of with household

waste.

Use proper landfill disposal or incineration in accordance

with local, state and federal regulations

**Uncleaned packaging:** 

**Recommendation:** Disposal must be made according to official regulations.

14. <u>Transport information</u>

Land transport ADR, RID: -

Maritime transport: IMDG-Code: -

Air transport ICAO-TI/IATA-DGR: -

15. Regulatory information

Classification according to EC-

guidelines:

The product is a medical device according to EC-

directive 93/42 EEC.

Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

# 16. Other information:

**Changes compared with the previous** Adaptation according to 1907/2006/EG, Article 31 **version:** 

The above information is based on our present day knowledge and relates solely to the safety requirements of the product. The data do not signify any warranty with regards to products properties. However users of the product should satisfy themselves that the information given is sufficient and correct for their specific circumstances of use.

Preparation by:	S. Zimmermann		
	Tel.: +49 521 8016800		