

NEW Renfert Polish hybrid materials

Diamond polishing paste

Diamond polishing paste specifically for high-gloss polishing of hybrid materials, for extra-oral use.



Renfert Polish hybrid materials

Diamond polishing paste specifically for high-gloss polishing of hybrid materials, for extra-oral use. Optimal polishing result thanks to diamond particles in the polishing paste that are specially matched to hybrid materials. Reduced plaque build-up and longer service life of the tooth restorations thanks to optimal surface finish.

Ordering information

Renfert Polish hybrid materials, 13 g

No. 510 3000

Further products of the Renfert Polish Series

Renfert Polish all-in-one

Universal polishing paste with a wide spectrum of diamond grain sizes, for extra-oral use. Quick, high-gloss finish on various materials such as lithium disilicate, zirconium dioxide, veneering ceramics and high performance polymers.

Renfert Polish ZrO₂

Extra-oral polishing paste for zirconium dioxide with specially coordinated abrasive properties due to extra coarse diamond grit.

Renfert Polish LiSi₂

Extra-oral polishing paste specially designed for lithium disilicate with a precise amount of coarse diamond grit.

Renfert Polish all-in-one, 18 g



No. 5100100

Renfert Polish ZrO₂, 18 g

No. 5101100

Renfert Polish LiSi₂, 18 g

No. 5102100






Renfert Polishing Pastes

Overview of Renfert Polishing Pastes with recommended Polishing Tools

Item		Recommended Polishing Tool				
Polishing Paste	with Dia- mond Grit	Material/ Range of Application	Processing/ Preliminary Polishing	Polishing	High Lustre Polishing	
Renfert Polish hybrid materials No. 510 3000	✓	Hybrid materials				
Renfert Polish all-in-one No. 5100100	✓	Lithium disilicate, zirconium dioxide, veneering ceramics, high performance polymers				
Renfert Polish ZrO ₂ No. 510 1100	✓	Zirconium dioxide				
Renfert Polish LiSi ₂ No. 510 2100	✓	Lithium disilicate				
Kohinoor L No. 5160001	✓	Veneering ceramic				
Brinell L No. 5190001	✓	Precious PFM alloys, ceramic				
Opal L No. 5200001		Veneering composites				
Saphir No. 515 0000		Precious metal alloys				
Saphir No. 515 0000		Cobalt-chrome alloys				
Universal polishing paste No. 513 1000		PMMA, resins				

*Included in Brinell L No. 5190001 scope of delivery

Polishers for Handpiece



Polishers for Lathe

