Safety	Data	Sheet
--------	------	-------

according to Regulation (EC) No 1907/2006

	GO 2011 speed	2012 0000
Revision date: 26.01.2021	Product code: 2012	20000 Page 1 of
SECTION 1: Identification of t	he substance/mixture and of the co	mpany/undertaking
<u>1.1. Product identifier</u> GO 2011 speed	2012 0000	
1.2. Relevant identified uses of t	he substance or mixture and uses advis	sed against
Use of the substance/mixture Plaster and Alginate Remo		
Uses advised against No information available.		
1.3. Details of the supplier of the		
Company name: Street: Place:	Renfert GmbH Untere Giesswiesen 2 D-78247 Hilzingen	
Telephone: e-mail: Contact person: e-mail: Internet:	+49 7731 8208-0 info@renfert.com Frau Andris silke.andris@renfert.com www.renfert.com	Telefax: +49 7731 8208-70
<u>1.4. Emergency telephone</u> number:	GB: 0044 151 951 3317 (Health and Safety Executive (HSE) Chemicals) USA: 800-222-1222 (National Poison Control number) CA: British Columbia: 800-567-8911 or 604-682-5050 Alberta: 1-800-332-1414 or 944-1414 Saskatchewan: 1-866-454-1212 Manitoba: 1-855-776-4766 Ontario: 1-800-268-9017 Quebec: 1-800-463-5060 New Brunswick: 911 New Foundland: 1-866-727-1110 Northwest Territories: 1-800-332-1414 Nova Scotia: 1-800-565-8161 AU: 13 11 26 ZA: 021 938 6084 (Office Hours) or 021 931 6129 (24 Hours)	
SECTION 2: Hazards identific	ation	

2.1. Classification of the substance or mixture

Regulation (EC) No. 1272/2008

Hazard categories: Skin corrosion/irritation: Skin Irrit. 2 Serious eye damage/eye irritation: Eye Irrit. 2 Hazard Statements: Causes skin irritation. Causes serious eye irritation.

2.2. Label elements

Renfert

Regulation (EC) No. 1272/2008

Signal word: Warning

Pictograms:



Hazard statements

H315	Causes skir
H319	Causes seri

Causes skin irritation. Causes serious eye irritation.

Precautionary statements

P280	Wear protective gloves and eye/face protection.
P337+P313	If eye irritation persists: Get medical advice/attention.

Renfert

Safety Data Sheet

according to Regulation (EC) No 1907/2006

GO 2011 speed 2012 0000

Revision date: 26.01.2021

Product code: 20120000

Page 2 of 8

2.3. Other hazards

The substances in the mixture do not meet the PBT/vPvB criteria according to REACH, annex XIII.

SECTION 3: Composition/information on ingredients

3.2. Mixtures

Chemical characterization

Complexing agent

Hazardous components

CAS No	Chemical name		Quantity	
	EC No	Index No	REACH No	
	GHS Classification	•		
65501-24-8	Tripotassium hydrogen ethylenedia	minetetraacetate		< 30 %
	241-543-5		01-2119904770-43	
	Acute Tox. 4; H332			
584-08-7	Potassium Carbonate			5 - 15 %
	209-529-3		01-2119532646-36	
	Skin Irrit. 2, Eye Irrit. 2, STOT SE 3	; H315 H319 H335		

Full text of H and EUH statements: see section 16.

Labelling for contents according to Regulation (EC) No 648/2004

15 % - < 30 % EDTA and salts thereof.

SECTION 4: First aid measures

4.1. Description of first aid measures

General information

Take off contaminated clothing and wash it before reuse.

After inhalation

Provide fresh air. When in doubt or if symptoms are observed, get medical advice.

After contact with skin

After contact with skin, wash immediately with plenty of water and soap. If skin irritation occurs: Get medical advice/attention.

After contact with eyes

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. In case of eye irritation consult an ophthalmologist.

After ingestion

Rinse mouth immediately and drink plenty of water. Medical treatment necessary.

4.2. Most important symptoms and effects, both acute and delayed

No information available.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media

The product itself does not burn. Co-ordinate fire-fighting measures to the fire surroundings.

Safety Data Sheet

according to Regulation (EC) No 1907/2006

GO 2011 speed 2012 0000

Revision date: 26.01.2021

Renfert

Product code: 20120000

Page 3 of 8

5.2. Special hazards arising from the substance or mixture

Thermal decomposition can lead to the escape of irritating gases and vapours.

5.3. Advice for firefighters

No special fire protection measures are necessary.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Use personal protection equipment.

6.2. Environmental precautions

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up

Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

6.4. Reference to other sections

Safe handling: see section 7 Personal protection equipment: see section 8 Disposal: see section 13

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Advice on safe handling

Avoid contact with eyes and skin. Wear personal protection equipment (refer to section 8). When using do not eat, drink or smoke.

Advice on protection against fire and explosion

No special fire protection measures are necessary.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

No special measures are necessary.

Further information on storage conditions

Do not store at temperatures below 0 °C

7.3. Specific end use(s)

Please refer to our internet website for more information: www.renfert.com

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Additional advice on limit values

To date, no national critical limit values exist.

8.2. Exposure controls

Appropriate engineering controls

See section 7. No additional measures necessary.

Protective and hygiene measures

Do not breathe vapour. Avoid contact with skin, eyes and clothes. Wash hands before breaks and after work. Keep away from food, drink and animal feedingstuffs.

Eye/face protection

Wear eye/face protection.

Renfert

Safety Data Sheet

according to Regulation (EC) No 1907/2006

GO 2011 speed 2012 0000

Revision date: 26.01.2021

Product code: 20120000

Page 4 of 8

Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. (DIN EN 374)

Suitable material: Butyl caoutchouc (butyl rubber), NBR (Nitrile rubber)

Thickness of the glove material: 0,4 mm

Breakthrough time: > 480 min

For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances.

Respiratory protection

Usually no personal respirative protection necessary. Provide adequate ventilation.

Environmental exposure controls

Do not allow to enter into surface water or drains.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state: Colour: Odour:	Liquid colourless odourless	
pH-Value (at 20 °C):	ododness	8
Changes in the physical state Melting point:		not determined
Initial boiling point and boiling range:		100 °C
Flash point:		not applicable
Flammability Solid: Gas:		not applicable not applicable
Explosive properties not explosive according to EU A.	14	
Lower explosion limits: Upper explosion limits:		not determined not determined
Ignition temperature:		not determined
Auto-ignition temperature Solid: Gas:		not applicable not applicable
Decomposition temperature: Oxidizing properties Not oxidising.		not determined
Vapour pressure: (at 20 °C)		23 hPa
Density (at 20 °C):		1,10 - 1,30 g/cm³
Water solubility:		completely miscible
Partition coefficient:		not determined
Viscosity / dynamic:		not determined
Viscosity / kinematic:		not determined

Safety Data Sheet

according to Regulation (EC) No 1907/2006

	GO 2011 speed	2012 0000	
Revision date: 26.01.2021	Product code	: 20120000	Page 5 of 8
Vapour density:		not determined	
Evaporation rate:		not determined	
9.2. Other information			
Solid content:		> 30 %	
No data available			

SECTION 10: Stability and reactivity

10.1. Reactivity

Renfert

No hazardous reaction when handled and stored according to provisions.

10.2. Chemical stability

The mixture is chemically stable under recommended conditions of storage, use and temperature.

10.3. Possibility of hazardous reactions

No known hazardous reactions.

10.4. Conditions to avoid

No special measures are necessary.

10.5. Incompatible materials

No information available.

10.6. Hazardous decomposition products

Carbon dioxide (CO2)

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity

Based on available data, the classification criteria are not met.

CAS No Chemical name

CAS NO	Chemical hame					
	Exposure route	Dose		Species	Source	Method
65501-24-8	Tripotassium hydrogen et	hylenediamir	netetraaceta	te		
	oral	LD50 mg/kg	2800	Rat	IUCLID	OECD 401
	inhalation vapour	ATE	11 mg/l			
	inhalation aerosol	ATE	1,5 mg/l			

Irritation and corrosivity

Causes skin irritation.

Causes serious eye irritation.

Sensitising effects

Based on available data, the classification criteria are not met.

Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

STOT-single exposure

Based on available data, the classification criteria are not met.

STOT-repeated exposure

Based on available data, the classification criteria are not met.

Aspiration hazard

Based on available data, the classification criteria are not met.

Renfert

Safety Data Sheet

according to Regulation (EC) No 1907/2006

GO 2011 speed 2012 0000

Revision date: 26.01.2021

Product code: 20120000

Page 6 of 8

Further information

The product has not been tested. The statement is derived from the properties of the single components.

SECTION 12: Ecological information

12.1. Toxicity

There are no data available on the mixture itself.

12.2. Persistence and degradability

There are no data available on the mixture itself.

12.3. Bioaccumulative potential

There are no data available on the mixture itself.

12.4. Mobility in soil

No information available.

12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to REACH, annex XIII.

12.6. Other adverse effects

No information available.

Further information

The product has not been tested. The statement is derived from the properties of the single components. Do not allow uncontrolled discharge of product into the environment.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Disposal recommendations

Dispose of waste according to applicable legislation.

The allocation of waste identity numbers/waste descriptions must be carried out according to the EEC, specific to the industry and process.

List of Wastes Code - residues/unused products

070701 WASTES FROM ORGANIC CHEMICAL PROCESSES; wastes from the MFSU of fine chemicals and chemical products not otherwise specified; aqueous washing liquids and mother liquors; hazardous waste

Contaminated packaging

Handle contaminated packages in the same way as the substance itself.

SECTION 14: Transport information

Land transport (ADR/RID)

Other applicable information (land transport)

No dangerous good in sense of these transport regulations.

Marine transport (IMDG)

Other applicable information (marine transport)

No dangerous good in sense of these transport regulations.

Air transport (ICAO-TI/IATA-DGR)

Other applicable information (air transport)

No dangerous good in sense of these transport regulations.

14.5. Environmental hazards

ENVIRONMENTALLY HAZARDOUS: no

14.6. Special precautions for user

not applicable

Safety Data Sheet

according to Regulation (EC) No 1907/2006

GO 2011 speed 2012 0000

Revision date: 26.01.2021

Renfert

Product code: 20120000

Page 7 of 8

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

EU regulatory information

Restrictions on use (REACH, annex XVII):

Entry 3: Flüssige Stoffe oder Gemische, die Kriterien für eine der folgenden in Anhang I der Verordnung (EG)

2004/42/EC (VOC):

not applicable

Additional information

Labelling for contents according to regulation (EC) No. 648/2004: 15 - < 30 % EDTA and salts thereof

National regulatory information

Water hazard class (D):

2 - obviously hazardous to water

Additional information

REGULATION (EC) No 1005/2009 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL on substances that deplete the ozone layer: not applicable REGULATION (EU) 2019/1021 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL on persistent organic pollutants : not applicable REGULATION (EU) No 649/2012 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL concerning the

export and import of hazardous chemicals: not applicable

SECTION 16: Other information

Changes

Abs. 11, 15, 16 * Data changed compared with the previous version

Abbreviations and acronyms

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

CAS: Chemical Abstracts Service (division of the American Chemical Society)

LD50: lethal dose, 50%

ADR: Accord européen relatif au transport international des marchandises Dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

RID: Règlement concernant le transport international ferroviaire de marchandises Dangereuses (Regulations

Concerning the International Transport of Dangerous Goods by Rail)

IMDG: International Maritime Code for Dangerous Goods

ICAO: International Civil Aviation Organization

MARPOL: International Convention for the Prevention of Marine Pollution from Ships

VOC: volatile organic compound(s)

Classification for mixtures and used evaluation method according to Regulation (EC) No. 1272/2008 [CLP]

Classification	Classification procedure
Skin Irrit. 2; H315	Calculation method
Eye Irrit. 2; H319	Calculation method

Relevant H and EUH statements (number and full text)

- H315 Causes skin irritation.
- H319 Causes serious eye irritation.
- H332 Harmful if inhaled.
- H335 May cause respiratory irritation.



according to Regulation (EC) No 1907/2006

GO 2011 speed 2012 0000		
	GO 2011 speed	2012 0000

Revision date: 26.01.2021

Renfert

Product code: 20120000

Page 8 of 8

Further Information

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material. Restricted to professional users.

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)