

**SECTION I - PRODUCT/ COMPANY IDENTIFICATION****PRODUCT NAME:** LUBECUBE

COMMON OR CHEMICAL NAME: Light-Cure Resin

PRODUCT (ITEM) NUMBER: 310108

**MANUFACTURED BY:** CENTRIX INC

EMERGENCY &amp; INFORMATION TELEPHONE: (800) 235-5862

Address: 770 RIVER ROAD

GENERAL INFORMATION: 203-929-5582

City, State, Zip: SHELTON, CT 06484

**SECTION II - HAZARDOUS INGREDIENTS/IDENTITY INFORMATION**

Product Contents: Light-Cure Bonding Agent Liquid

Hazardous Ingredients: None listed.

**SECTION III - HAZARDS IDENTIFICATION**

For potential health hazards: SEE Section 11 below.

No other significant hazards to man or the environment are reasonably foreseen.

**SECTION IV - FIRST AID MEASURES**

Inhalation: No harm anticipated.

Skin Contact: Wash with soap and water.

Eye Contact: Flush with copious amounts of water.

Ingestion: Induce vomiting if large amount of material

**SECTION V - FIRE-FIGHTING MEASURES**

Suitable Extinguishing Media: Water contact acceptable.

Extinguishing Media which should not be used: None.

Special hazards from combustion products: None specific.

Special protective equipment for fire-fighters: None.

**SECTION VI - ACCIDENTAL RELEASE MEASURES**

Personal Precautions: Avoid contact of skin or eyes.

Wash hands with soap and water if contact occurs.

Environmental Precautions: Dispose in accordance with all federal, state, and local regulations.

Methods for Cleaning Up: Wipe up with soap, water, and paper towels.

**SECTION VII - HANDLING AND STORAGE**

Handling Precautions: Not applicable.

Safe Storage Conditions: Store at temperatures not exceeding 75°F (24°C).

Other Restrictions or Recommendations: Avoid exposure to moisture and light.

**SECTION VIII - EXPOSURE CONTROLS / PERSONAL PROTECTION**

Respiratory Protection &amp; Ventilation: Not required.

Protective Clothing or Equipment: Not required under normal handling.

Unusual Work/Hygienic Practices: Avoid prolonged contact with skin; wash with soap and water.

**SECTION IX - PHYSICAL AND CHEMICAL PROPERTIES**

Appearance: Clear Liquid.

Odor: Slight acrylic odor.

pH: Not applicable.

Boiling Point/Range: Not applicable.

Melting Point/Range: Not applicable.

Flash Point: Not applicable.

Auto-Ignition Temperature: Not applicable.

Flammability: No.

Explosive Limits in Air: Not applicable

**(See next column)**

Oxidizing Properties: None.

Vapor pressure (mm Hg): Negligible at room temperatures.

Vapor Density (Air=1): &gt;1

Specific Gravity (H2O=1): ~1

Solubility in Water: Insoluble.

Reactivity in Water: Not reactive.

Other Data: None.

## **SECTION X - STABILITY AND REACTIVITY**

Stability: Generally stable.

Conditions of Instability: Avoid exposure to light.

Incompatibility with Various Substances: None

Hazardous Decomposition Products: None

Hazardous Polymerization: May occur upon exposure to light where large quantities of material are involved.

Conditions to Avoid: Avoid exposure to light and temperatures above 80°F.

## **SECTION XI - TOXICOLOGICAL INFORMATION**

Potential Acute Health Effects: Possibility of skin irritation or sensitization if not washed off promptly.

Potential Chronic Health Effects: Possibility of skin sensitization or irritation upon prolonged or repeated exposure.

Signs and Symptoms of Exposure: Skin/eye irritation.

Medical Conditions Generally Aggravated by Exposure: None known.

Carcinogenic or Potential Carcinogenic Effects: No ingredient listed as a carcinogen.

Other Potential Health Effects: None known.

## **SECTION XII - ECOLOGICAL INFORMATION**

No ecological damage or impact anticipated when handled as indicated.

## **SECTION XIII - DISPOSAL CONSIDERATIONS**

Waste Disposal: Do not dispose into sink.

Dispose of in accordance with all federal, state, and local regulations.

## **SECTION XIV - TRANSPORT INFORMATION**

UN Classification: Not applicable.

Identification: Not applicable.

Special Provisions for Transport: Avoid elevated temperatures and high intensity light.

## **SECTION XV - REGULATORY INFORMATION**

U.S. FDA Dental Device: Class II product

EC Medical Device Directive: Class IIa product

## **SECTION XVI - OTHER INFORMATION**

Some of the information presented and conclusions drawn herein are from sources other than direct test data on the product itself. The information in this MSDS was obtained from sources that we believe are reliable. However, the information is provided without any warranty, express or implied, regarding its correctness. The conditions or method of handling, storage, use and disposal of the product are beyond our control and may be beyond our knowledge. For this and other reasons, we do not assume responsibility and expressly disclaim liability for loss, damage, or expense arising out of or in any way connected with the handling, storage, use or disposal of the product. If the product is used as a component in another product or used in a way other than recommended by the Company, this MSDS information may not be applicable.

**centrix**<sup>®</sup>

*Making Dentistry Easier.<sup>SM</sup>*

### **CENTRIX INCORPORATED**

770 River Road, Shelton, CT 06484 USA  
(800) 235-5862 Tel: (203) 929-5582  
Fax: (203) 929-6804

**CE**  
0510

<b>EC</b>	<b>REP</b>
-----------	------------

### **CENTRIX EUROPE**

Hauptstrasse 64 B, 65719 Hofheim, Germany

©Centrix, Inc. 2008 Date Revised 10/19/2010