

Safety data sheet according to 1907/2006/EC, Article 31

Printing date 04.08.2023

Version number 4

Revision date: 04.08.2023

SECTION 1: Identification of the substance/mixture and of the company/undertaking

· 1.1 Product identifier

[·] Trade name: PANAVIA V5 ; Try-in Paste

- · Article number: 214-Ti
- **Registration number** Not applicable
- 1.2 Relevant identified uses of the substance or mixture and uses advised against
- No further relevant information available.
- · Application of the substance / the mixture Dental colour matching test material

· 1.3 Details of the supplier of the safety data sheet

- Supplier : Kuraray Europe GmbH Philipp-Reis-Str. 4, 65795 Hattersheim am Main Germany Phone: +49 (0)69 305 35 835 E-mail: dental.eu@kuraray.com
- · Further information obtainable from: Supplier
- 1.4 Emergency telephone number:
- 3E Company ; (UK) 0 800 680 0425
 - (GB) +44 20 35147487 (AU) +61 2 80363166
 - Kuraray Access Code: 334674

SECTION 2: Hazards identification

· 2.1 Classification of the substance or mixture

· Classification according to Regulation (EC) No 1272/2008

Carc. 2 H351 Suspected of causing cancer. Route of exposure: Inhalation.

· Additional information

Medical devices as defined in Regulation (EU) 2017/745 and which are invasive or used in direct physical contact with the human body, are exempted from the provisions of Regulation (EC) No 1272/2008 (GB CLP/ GHS) usually if they are in the finished state and intended for the final user.

· 2.2 Label elements

· Labelling according to Regulation (EC) No 1272/2008

- The product is classified and labelled according to the GB CLP regulation.
- Hazard pictograms



· Signal word Warning

· Hazard-determining components of labelling:

- titanium dioxide
- · Hazard statements
- H351 Suspected of causing cancer. Route of exposure: Inhalation.

· Precautionary statements

- P201 Obtain special instructions before use.
- P202 Do not handle until all safety precautions have been read and understood.
- P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.
- P308+P313 IF exposed or concerned: Get medical advice/attention.
- P405 Store locked up.
- P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

(Contd. on page 2)

GB



Safety data sheet according to 1907/2006/EC, Article 31

Printing date 04.08.2023

Version number 4

Trade name: PANAVIA V5 ; Try-in Paste

· Labelling of packages where the contents do not exceed 125 ml

· Precautionary statements

EUH212: Warning! Hazardous respirable dust may be formed when used. Do not breathe dust.

· 2.3 Other hazards

· Results of PBT and vPvB assessment

- **PBT:** Not applicable.
- **vPvB:** Not applicable.

SECTION 3: Composition/information on ingredients

• Description: Mixture of substances listed below with nonhazardous additions.

· Dangerous components:				
CAS: 56-81-5 EINECS: 200-289-5	Glycerol	substance with a Community workplace exposure limit	75-95%	
CAS: 60676-86-0 EINECS: 262-373-8		substance with a Community workplace exposure limit	<10%	
CAS: 13463-67-7 EINECS: 236-675-5	titanium dioxide	& Carc. 2, H351	<5%	

Other ingredients:

Silanated amorphous silica

- Amorphous silica
- Pigments
- Additional information:

Silicon oxide; Nanoform, Amorphous silica and Silanated amorphous silica, particle size 1 - 100 nm

SECTION 4: First aid measures

· 4.1 Description of first aid measures

- After inhalation: Supply fresh air and to be sure call for a doctor.
- After skin contact: Generally the product does not irritate the skin.
- · After eye contact: Rinse opened eye for several minutes under running water.
- After swallowing: If symptoms persist consult doctor.
- · 4.2 Most important symptoms and effects, both acute and delayed

No further relevant information available.

· 4.3 Indication of any immediate medical attention and special treatment needed

No further relevant information available.

SECTION 5: Firefighting measures

- · 5.1 Extinguishing media
- Suitable extinguishing agents: Use fire extinguishing methods suitable to surrounding conditions.
- 5.2 Special hazards arising from the substance or mixture No further relevant information available.
- 5.3 Advice for firefighters
- · Protective equipment: No special measures required.

SECTION 6: Accidental release measures

- 6.1 Personal precautions, protective equipment and emergency procedures See Section 8.
- \cdot 6.2 Environmental precautions: Do not allow to enter sewers/ surface or ground water.
- **6.3** Methods and material for containment and cleaning up:

Pick up mechanically. Ensure adequate ventilation.

(Contd. on page 3)

(Contd. of page 1)

^{· 3.2} Mixtures



Printing date 04.08.2023

Safety data sheet according to 1907/2006/EC, Article 31

Version number 4

Revision date: 04.08.2023

Trade name: PANAVIA V5 ; Try-in Paste

(Contd. of page 2)

• 6.4 Reference to other sections No dangerous substances are released.

SECTION 7: Handling and storage

· 7.1 Precautions for safe handling No special measures required.

· Information about fire - and explosion protection: No special measures required.

· 7.2 Conditions for safe storage, including any incompatibilities

· Storage:

• Requirements to be met by storerooms and receptacles: Store in a cool location (2-25°C).

· Information about storage in one common storage facility: Not required.

- Further information about storage conditions: Store in cool, dry conditions in well sealed receptacles.
- Keep container tightly sealed.

Protect from heat and direct sunlight.

• 7.3 Specific end use(s) No further relevant information available.

SECTION 8: Exposure controls/personal protection

• 8.1 Control parameters

· 8.1 Control parameters			
· Ingredients with limit values that require monitoring at the workplace:			
56-81-5 Glycerol			
WEL	Long-term value: 10 mg/m ³		
60676-86-0 Silica, fused			
WEL (Great Britain)	Long-term value: 0.08 mg/m ³		
	fused respirable dust		
13463-67-7 titanium dioxide			
WEL (Great Britain)	Long-term value: 10* 4** mg/m ³		
	*total inhalable **respirable		
· Additional information: The lists valid during the making were used as basis.			

· 8.2 Exposure controls

• Appropriate engineering controls No further data; see section 7.

· Personal protective equipment:

· General protective and hygienic measures:

The usual precautionary measures are to be adhered to when handling chemicals.

· Respiratory protection: Suitable respiratory protective device recommended.

· Hand protection

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material

The exact break trough time has to be found out by the manufacturer of the protective gloves and has to be observed.

(Contd. on page 4)

GB -

(Contd. of page 3)



Printing date 04.08.2023

*

Safety data sheet according to 1907/2006/EC, Article 31

Version number 4

Revision date: 04.08.2023

Trade name: PANAVIA V5 ; Try-in Paste

 \cdot Eye/face protection Goggles recommended during refilling

SECTION 9: Physical and chemical properties

• 9.1 Information on basic physical and chemical pr	-		
· Physical state	Solid		
· Colour:	Coloured		
· Odour:	Odourless		
 Melting point/freezing point: 	18.2 °C (glycerol)		
· Boiling point or initial boiling point and boiling			
range	290.0 °C (glycerol)		
· Flammability	Not determined.		
· Lower and upper explosion limit			
· Lower:	0.9 Vol % (glycerol)		
· Upper:	Not determined		
· Flash point:	160 °C (glycerol)		
· Auto-ignition temperature:	400.0 °C (glycerol)		
· Decomposition temperature:	Not determined		
	Not determined		
· pH · Salukility			
· Solubility	Dorthy soluble		
· water:	Partly soluble.		
· Partition coefficient n-octanol/water (log value)	Not determined		
Vapour pressure:	0.1 hPa (20°C, glycerol)		
Density and/or relative density			
Density:	1.2 g/cm^3		
· Relative density	Not determined		
· Relative vapour density	Not determined		
 Particle characteristics 			
See section 3.			
• 9.2 Other information			
· Important information on protection of health an	d		
environment, and on safety.			
· Explosive properties:	Product does not present an explosion hazard.		
· Change in condition	r roudet does not present un expresion nuzura.		
· Softening point/range			
· Oxidising properties	Not determined		
· Oxidising gases	Void		
	Void		
· Gases under pressure	Void		
· Flammable liquids · Flammable solids	Void		
• Self-reactive substances and mixtures	Void Void		
· Pyrophoric liquids	Void		
· Pyrophoric solids	Void		
Self-heating substances and mixtures	Void		
Substances and mixtures, which emit flammable			
gases in contact with water	Void		
· Oxidising liquids	Void		
· Oxidising solids	Void		
· Organic peroxides	Void		
· Corrosive to metals	Void		
· Desensitised explosives	Void		

(Contd. on page 5)



Printing date 04.08.2023

Safety data sheet according to 1907/2006/EC, Article 31

Version number 4

Revision date: 04.08.2023

Trade name: PANAVIA V5 ; Try-in Paste

(Contd. of page 4)

SECTION 10: Stability and reactivity

- · 10.1 Reactivity No further relevant information available.
- · 10.2 Chemical stability
- Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- 10.3 Possibility of hazardous reactions No dangerous reactions known.
- 10.4 Conditions to avoid No further relevant information available.
- 10.5 Incompatible materials: No further relevant information available.
- 10.6 Hazardous decomposition products: Carbon monoxide and carbon dioxide

SECTION 11: Toxicological information

- · 11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008
- Acute toxicity Based on available data, the classification criteria are not met.
- · Skin corrosion/irritation Based on available data, the classification criteria are not met.
- · Serious eye damage/irritation Based on available data, the classification criteria are not met.
- Respiratory or skin sensitisation Based on available data, the classification criteria are not met.
- · Germ cell mutagenicity Based on available data, the classification criteria are not met.
- · Carcinogenicity
- Suspected of causing cancer. Route of exposure: Inhalation.
- · Reproductive toxicity Based on available data, the classification criteria are not met.
- STOT-single exposure Based on available data, the classification criteria are not met.
- · STOT-repeated exposure Based on available data, the classification criteria are not met.
- · Aspiration hazard Based on available data, the classification criteria are not met.
- · 11.2 Information on other hazards

· Endocrine disrupting properties

None of the ingredients is listed.

SECTION 12: Ecological information

- · 12.1 Toxicity
- · Aquatic toxicity: No further relevant information available.
- 12.2 Persistence and degradability No further relevant information available.
- 12.3 Bioaccumulative potential No further relevant information available.
- 12.4 Mobility in soil No further relevant information available.
- · 12.5 Results of PBT and vPvB assessment
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- · 12.6 Endocrine disrupting properties
- The product does not contain substances with endocrine disrupting properties.
- · 12.7 Other adverse effects
- · Additional ecological information:
- · General notes:

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

SECTION 13: Disposal considerations

- · 13.1 Waste treatment methods
- · Recommendation Smaller quantities can be disposed of with household waste.

(Contd. on page 6)



Safety data sheet according to 1907/2006/EC, Article 31

Printing date 04.08.2023

Version number 4

Revision date: 04.08.2023

Trade name: PANAVIA V5 ; Try-in Paste

(Contd. of page 5)

- · Uncleaned packaging:
- Recommendation: Disposal must be made according to official regulations.

SECTION 14: Transport information

SECTION 14. Transport mitormation		
· 14.1 UN number or ID number · ADR, ADN, IMDG, IATA	Void	
· 14.2 UN proper shipping name · ADR, ADN, IMDG, IATA	Void	
· 14.3 Transport hazard class(es)		
· ADR, ADN, IMDG, IATA · Class	Void	
· 14.4 Packing group · ADR, ADN, IMDG, IATA	Void	
· 14.5 Environmental hazards: · Marine pollutant:	No	
· 14.6 Special precautions for user	Not applicable.	
• 14.7 Maritime transport in bulk according to IMO instruments Not applicable.		
· UN "Model Regulation":	Void	

SECTION 15: Regulatory information

· 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

- Directive 94/33/EC: Care should be taken when people under 18 handle this product.
- Directive 94/85/EEC: No further relevant information available.
- · Named dangerous substances ANNEX I None of the ingredient is listed.

· Regulation (EU) No 649/2012

· Regulation (EC) No 273/2004 on drug precursors

None of the ingredients is listed.

• Regulation (EC) No 111/2005 laying down rules for the monitoring of trade between the Community and third countries in drug precursors

None of the ingredients is listed.

15.2 Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Relevant phrases

H351 Suspected of causing cancer.

• Abbreviations and acronyms:

ELINCS: European List of Notified Chemical Substances

Carc. 2: Carcinogenicity – Category 2

 \cdot * Data compared to the previous version altered.