



SAFETY DATA SHEET (SDS)

In accordance with Hazard Communication Standard (HCS), 29CFR 1910.1200(g), 2012 revision

SECTION I - PRODUCT / COMPANY IDENTIFICATION

Product Trade Name: NoMIX® Temporary Cement.
Description: Moisture activated temporary cement for crowns and bridges
Part (Item) Numbers: 290040, 310157, 310520
Company Name: Centrix, Inc.
Address: 770 River Road
City, State, Zip: Shelton CT, 06484
Emergency Telephone Number: (203) 929-5582
Information Telephone Number: (800) 235-5862
Date Revised: 09/2014

SECTION II - HAZARD IDENTIFICATION

Designation of Risk: Xi, irritant  
Special Details of Risk (R-categories): 36-38
Detrimental effects/symptoms: Irritant in case of eye, skin and inhalation contact.
May cause sensitization by skin contact.

SECTION III - COMPOSITION / INFORMATION ON INGREDIENTS

Component	%	CAS#	EINECS #
D.M.T.G (Glycol)	<20%	143-24-8	205-594-7
Ethyl Methacrylate Polymer	<4%		
Calcium Sulfate	<31%	10034-76-1	231-900-3
Zinc Oxide	<26%	1314-13-2	
Triacetin	<3%	102-76-1	
Potassium Sulfate	<21%	7778-80-5	231-915-6

Exact % content withheld due to trade secret restrictions.

SECTION IV - FIRST AID MEASURES

Skin: Wash with soap and water.
Ingestion: N/A Inhalation: N/A
Eyes: Flush with copious amounts of water.
If irritation persists seek medical attention

SECTION V - FIRE-FIGHTING MEASURES

Burning produces irritating, toxic and obnoxious fumes.
Suitable extinguishing media: Use extinguisher suitable to cause of fire and its surroundings.
Forbidden extinguishing media: not available
Special firefighting measures: not available
Unusual fire and explosion hazards: not available
Special protection equipment: Wear suitable respiratory equipment when necessary.

SECTION VI - ACCIDENTAL RELEASE MEASURES

Transfer bulk into suitable containers for disposal. Wipe residue with damp cloth and scrape up into suitable container for disposal.
Personal precautions: Ensure adequate ventilation of the working area.
Environmental precautions: Do not allow product to enter drains. Prevent further spillage if safe.
Reclaiming methods: Sweep up. Transfer to suitable, labelled containers for disposal.
Clean spillage area thoroughly with plenty of water.

SECTION VII - HANDLING AND STORAGE

Handling: No handling problems are likely to be encountered provided normal safe practices of hygiene are applied.
Storage: Store in a cool dry place.
Containers: Keep containers tightly closed. Store in correctly labelled containers
Incompatibilities: unknown
Other precautions: no further information available

SECTION VIII - EXPOSURE CONTROLS / PERSONAL PROTECTION

General measures: Wash hands thoroughly after contact with this product.
Limits: no data is available
Respiratory protection: Wear suitable respiratory equipment when necessary.
Eye protection: Goggles / eyeshields
Skin protection: Overalls, rubber gloves
Hands protection: Overalls, rubber gloves
Other equipment: no further information available

SECTION IX - PHYSICAL AND CHEMICAL PROPERTIES

Appearance: White Stiff Paste
Odor: Odorless
Boiling point: N/A
Flash point N/A
Other Data: Not available or none.
Explosion limits: N/A
Water solubility: Insoluble
Evaporation rate: N/A

SECTION X - STABILITY AND REACTIVITY.

Stability: Stable, non-reactive
Conditions to avoid: N/A
Hazardous decomposition products: N/A
Materials to avoid: N/A
Hazardous polymerization products: N/A

SECTION XI - TOXICOLOGICAL INFORMATION.

No data is available

SECTION XII - ECOLOGICAL INFORMATION.

No data is available

SECTION XIII - DISPOSAL CONSIDERATIONS.

No data is available

SECTION XIV - TRANSPORT INFORMATION.

No data is available

SECTION XV - REGULATORY INFORMATION

No data is available

SECTION XVI - OTHER INFORMATION

Handling of product exclusively by dental specialists according to instructions for use.

CAUTION: The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. Any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user's responsibility.

centrix®

Making Dentistry Easier. SM



CENTRIX INCORPORATED

770 River Road, Shelton, CT 06484 USA
(203) 929-5582 • (800) 235-5862 • FAX (203) 929-6804
www.centrixdental.com



EC REP Dental-Contact Vertriebs OHG, Im Spiet 101, DE 26506 Norden

EMERGO AUSTRALIA

Level 20, Tower II, Darling Park, 201 Sussex Street
Sydney, NSW 2000, Australia