

Safety Data Sheet

according to UK REACH Regulation

GB600 STERILIT Power Systems Oil spray

Revision date: 20.02.2024

Product code: 00250-0019

Page 1 of 10

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

GB600 STERILIT Power Systems Oil spray

UFI: TTSA-H6KK-X70D-22SP

1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/mixture

Preservative agent

1.3. Details of the supplier of the safety data sheet

Company name: AESCULAP AG
Street: Am AESCULAP-Platz
Place: D-78532 Tuttlingen/Donau
Telephone: +49 (0) 74 61 95-0
Responsible Department: Responsible for the safety data sheet: sds@gbk-ingelheim.de
Telefax: +49 (0) 74 61 95-2600

1.4. Emergency telephone number:

INTERNATIONAL: +49 - (0) 6132 - 84463, GBK GmbH (24h - 7d/w - 365d/a)
In England and Wales: NHS 111 In Scotland: NHS 24 - dial 111

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

GB CLP Regulation

Aerosol 1; H222-H229
Aquatic Chronic 3; H412

Full text of hazard statements: see SECTION 16.

2.2. Label elements

GB CLP Regulation

Hazard components for labelling

Pentane

Signal word: Danger

Pictograms:



Hazard statements

H222 Extremely flammable aerosol.
H229 Pressurised container: May burst if heated.
H412 Harmful to aquatic life with long lasting effects.

Precautionary statements

P102 Keep out of reach of children.
P210 Keep away from heat. No Smoking.
P211 Do not spray on an open flame or other ignition source.
P251 Do not pierce or burn, even after use.
P410+P412 Protect from sunlight. Do not expose to temperatures exceeding 50 °C/122 °F.

2.3. Other hazards

Safety Data Sheet

according to UK REACH Regulation

GB600 STERILIT Power Systems Oil spray

Revision date: 20.02.2024

Product code: 00250-0019

Page 2 of 10

According to Regulation (EC) No 1907/2006 (REACH) none of the substances, contained in this product are a PBT / vPvB substance.

The mixture does not contain substance(s) included in the list established in accordance with Article 59(1) of REACH for having endocrine disrupting properties, or is not identified as having endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605

In use, may form flammable/explosive vapour-air mixture.

SECTION 3: Composition/information on ingredients**3.2. Mixtures****Chemical characterization**

Aerosol with preparation of special oil and pentane

Volume: 300 ml

Propellant: propane/butan < 150 g

Hazardous components

| CAS No | Chemical name | | | Quantity |
|-----------|---|--------------|------------------|------------|
| | EC No | Index No | REACH No | |
| | Classification (GB CLP Regulation) | | | |
| 106-97-8 | Butane | | | 60 - 100 % |
| | 203-448-7 | 601-004-00-0 | 01-2119474691-32 | |
| | Flam. Gas 1, Press. Gas (Comp.); H220 H280 | | | |
| 74-98-6 | Propane | | | 15 - 45 % |
| | 200-827-9 | 601-003-00-5 | 01-2119486944-21 | |
| | Flam. Gas 1, Press. Gas (Comp.); H220 H280 | | | |
| 109-66-0 | Pentane | | | < 15 % |
| | 203-692-4 | 601-006-00-1 | 01-2119459286-30 | |
| | Flam. Liq. 1, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2; H224 H336 H304 H411 EUH066 | | | |
| 8042-47-5 | White mineral oil (petroleum) | | | < 10 % |
| | 232-455-8 | | 01-2119487078-27 | |
| | Asp. Tox. 1; H304 | | | |

Full text of H and EUH statements: see section 16.

SECTION 4: First aid measures**4.1. Description of first aid measures****After inhalation**

Ensure of fresh air.

In the event of symptoms refer for medical treatment.

After contact with skin

Wash with water and soap as a precaution.

Consult a doctor if skin irritation persists.

After contact with eyes

Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes.

If eye irritation persists, consult a specialist.

After ingestion

Do not induce vomiting.

Summon a doctor immediately.

Safety Data Sheet

according to UK REACH Regulation

GB600 STERILIT Power Systems Oil spray

Revision date: 20.02.2024

Product code: 00250-0019

Page 3 of 10

Induce vomiting only upon the advice of a physician.

4.2. Most important symptoms and effects, both acute and delayed

Inhalation of high vapour concentrations may cause symptoms like headache, dizziness, tiredness, nausea and vomiting.

Evaporative cooling effect during use.

Contact with eyes, skin or oral tissues may cause irritation.

Attention. Beware, danger of aspiration.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptoms.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media

Foam, carbon dioxide (CO₂), dry chemical, water-spray.

Unsuitable extinguishing media

Full water jet

5.2. Special hazards arising from the substance or mixture

Fire may produce:

carbon monoxide and carbon dioxide

5.3. Advice for firefighters

Use breathing apparatus with independent air supply.

Protective suit.

Additional information

Cool containers at risk with water spray jet. Heating will cause pressure rise with risk of bursting.

Vapours are heavier than air and spread along ground.

The vapour/air mixture is explosive, even in empty, uncleaned receptacles.

Fire residues and contaminated firefighting water must be disposed of in accordance with the local regulations.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General advice

Ensure adequate ventilation.

Remove all sources of ignition.

Keep away noninvolved persons.

For non-emergency personnel

Do not breath vapours or spray mist.

Avoid contact with skin, eyes and clothing.

For emergency responders

Use personal protective clothing.

Use breathing apparatus if exposed to vapours/aerosol.

6.2. Environmental precautions

Do not discharge into the drains/surface waters/ground water.

Do not discharge into the subsoil/soil.

6.3. Methods and material for containment and cleaning up

For containment

Prevent spread over a wide area (e.g. by containment or oil barriers).

For cleaning up

Soak up with inert absorbent material (e.g. sand, silica gel, acid binder, universal binder).

Shovel into suitable container for disposal.

Safety Data Sheet

according to UK REACH Regulation

GB600 STERILIT Power Systems Oil spray

Revision date: 20.02.2024

Product code: 00250-0019

Page 4 of 10

6.4. Reference to other sections

Observe protective instructions (see Sections 7 and 8).
Information for disposal see section 13.

SECTION 7: Handling and storage**7.1. Precautions for safe handling****Advice on safe handling**

Ensure adequate ventilation.
Avoid contact with skin, eyes and clothing.
Do not breathe aerosol

Advice on protection against fire and explosion

Pressurized container: protect from sunlight and do not expose to temperatures exceeding 50 °C. Do not pierce or burn, even after use.
Do not spray on a naked flame or any incandescent material. Keep away from sources of ignition - No smoking. Keep out of the reach of children. Without adequate ventilation formation of explosive mixtures may be possible.

Advice on general occupational hygiene

Do not breathe aerosol
Wash hands before breaks and at the end of workday.
When using do not eat, drink or smoke.

7.2. Conditions for safe storage, including any incompatibilities**Requirements for storage rooms and vessels**

Store at ambient temperature in a well-ventilated place. Keep container tightly closed.

Hints on joint storage

Incompatible with oxidizing agents.

Further information on storage conditions

Keep away from food, drink and animal feeding stuffs.

7.3. Specific end use(s)

Preservative agent

SECTION 8: Exposure controls/personal protection**8.1. Control parameters****Exposure limits (EH40)**

| CAS No | Substance | ppm | mg/m ³ | fibres/ml | Category | Origin |
|----------|-----------|-----|-------------------|-----------|---------------|--------|
| 106-97-8 | Butane | 600 | 1450 | | TWA (8 h) | WEL |
| | | 750 | 1810 | | STEL (15 min) | WEL |
| 109-66-0 | Pentane | 600 | 1800 | | TWA (8 h) | WEL |

DNEL/DMEL values

| CAS No | Substance | Exposure route | Effect | Value |
|--------------------------|-------------------------------|----------------|----------|-----------------------|
| 8042-47-5 | White mineral oil (petroleum) | | | |
| Consumer DNEL, long-term | | inhalation | systemic | 160 mg/m ³ |
| Consumer DNEL, long-term | | dermal | systemic | 220 mg/kg bw/day |

8.2. Exposure controls

Safety Data Sheet

according to UK REACH Regulation

GB600 STERILIT Power Systems Oil spray

Revision date: 20.02.2024

Product code: 00250-0019

Page 5 of 10

Appropriate engineering controls

Ensure adequate ventilation, especially in confined areas.

Individual protection measures, such as personal protective equipment**Eye/face protection**

If used properly, no need to wear eye protection.

Avoid contact with the eyes.

Hand protection

If used properly, protective gloves are normally not required.

Avoid contact with the skin.

Respiratory protection

No personal respiratory protective equipment normally required.

In case of insufficient ventilation wear suitable respiratory equipment (gas filter type AX) (EN 14387).

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties**

| | |
|-----------------|----------------|
| Physical state: | Aerosol |
| Colour: | Colourless |
| Odour: | Characteristic |

Test method**Changes in the physical state**

| | |
|---|---------------|
| Melting point/freezing point: | < - 138 °C *) |
| Boiling point or initial boiling point and boiling range: | < 0 °C *) |
| Sublimation point: | n.a. |
| Softening point: | n.a. |
| Flash point: | - 80 °C *) |

Flammability

| | |
|---------------|------|
| Solid/liquid: | n.a. |
|---------------|------|

Explosive properties

In use, may form flammable/explosive vapour-air mixture.

| | |
|----------------------------|-----------------|
| Lower explosion limits: | 1,5 vol. % *) |
| Upper explosion limits: | 10,9 vol. % *) |
| Auto-ignition temperature: | 365 - 470 °C *) |

Self-ignition temperature

| | |
|--------|------|
| Solid: | n.a. |
| Gas: | n.a. |

| | |
|----------------------------|------|
| Decomposition temperature: | n.d. |
|----------------------------|------|

| | |
|-----------|----------------------------|
| pH-Value: | Not detectable by analysis |
|-----------|----------------------------|

| | |
|----------------------|------|
| Viscosity / dynamic: | n.a. |
|----------------------|------|

| | |
|------------------------|------|
| Viscosity / kinematic: | n.a. |
|------------------------|------|

| | |
|------------|------|
| Flow time: | n.a. |
|------------|------|

| | |
|---------------------------------|------------|
| Water solubility: (at 20 °C) | Immiscible |
|---------------------------------|------------|

Solubility in other solvents

n.d.

| | |
|--|------|
| Partition coefficient n-octanol/water: | n.d. |
|--|------|

Safety Data Sheet

according to UK REACH Regulation

GB600 STERILIT Power Systems Oil spray

Revision date: 20.02.2024

Product code: 00250-0019

Page 6 of 10

| | |
|--------------------------------|------------------------|
| Vapour pressure: (at 20 °C) | 2200 - 8400 hPa *) |
| Density (at 20 °C): | 0,58 g/cm ³ |
| Bulk density: | n.a. |
| Relative vapour density: | n.d. |

9.2. Other information**Information with regard to physical hazard classes**

| | |
|--|-----------------------|
| Sustaining combustion: | Sustaining combustion |
| Oxidizing properties Not oxidising. | |

Other safety characteristics

| | |
|-------------------|------|
| Evaporation rate: | n.d. |
|-------------------|------|

Further Information

*) Propellant

SECTION 10: Stability and reactivity**10.1. Reactivity**

No decomposition if stored and applied as directed.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

Reactions with oxidizing agents.

10.4. Conditions to avoid

Fire or intense heat may cause violent rupture of packages.

10.5. Incompatible materials

Strong oxidizing agents

10.6. Hazardous decomposition products

No hazardous decomposition products known.

Fire may produce:
carbon monoxide and carbon dioxide**SECTION 11: Toxicological information****11.1. Information on hazard classes as defined in GB CLP Regulation****Acute toxicity**Based on available data, the classification criteria are not met.
No toxicological data available.**ATEmix calculated**

ATE (oral) > 2000 mg/kg; ATE (dermal) > 2000 mg/kg; ATE (inhalation vapour) > 20 mg/l; ATE (inhalation dust/mist) > 5 mg/l

Irritation and corrosivity

Based on available data, the classification criteria are not met.

Sensitising effects

Based on available data, the classification criteria are not met.

Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

STOT-single exposure

Based on available data, the classification criteria are not met.

Safety Data Sheet

according to UK REACH Regulation

GB600 STERILIT Power Systems Oil spray

Revision date: 20.02.2024

Product code: 00250-0019

Page 7 of 10

STOT-repeated exposure

Based on available data, the classification criteria are not met.

Aspiration hazard

Based on available data, the classification criteria are not met.

Additional information on tests

Classification in compliance with the assessment procedure specified in the Regulation (EC) no 1272/2008.

11.2. Information on other hazards**Endocrine disrupting properties**

No data available

Other information

Inhalation of high vapour concentrations may cause symptoms like headache, dizziness, tiredness, nausea and vomiting.

Evaporative cooling effect during use.

Contact with eyes, skin or oral tissues may cause irritation.

SECTION 12: Ecological information**12.1. Toxicity**

Ecological data are not available.

Harmful to aquatic life with long lasting effects.

12.2. Persistence and degradability

No data available

| CAS No | Chemical name | | | |
|-----------|-------------------------------|--------|----|--------|
| | Method | Value | d | Source |
| | Evaluation | | | |
| 8042-47-5 | White mineral oil (petroleum) | | | |
| | OECD 301 F | 31,3 % | 28 | |

12.3. Bioaccumulative potential

No data available

12.4. Mobility in soil

No data available

12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to UK REACH.

12.6. Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

12.7. Other adverse effects

Hazardous water pollutant.

Further information

Do not flush into surface water or sanitary sewer system.

Ecological injuries are not known or expected under normal use.

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Disposal recommendations**

Can be incinerated, when in compliance with local regulations.

Where possible recycling is preferred to disposal.

List of Wastes Code - residues/unused products

Safety Data Sheet

according to UK REACH Regulation

GB600 STERILIT Power Systems Oil spray

Revision date: 20.02.2024

Product code: 00250-0019

Page 8 of 10

150111 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); metallic packaging containing a hazardous solid porous matrix (for example asbestos), including empty pressure containers; hazardous waste

Contaminated packaging

Offer empty spray cans to an established disposal company.

SECTION 14: Transport information

Land transport (ADR/RID)

14.1. UN number or ID number: UN 1950
14.2. UN proper shipping name: AEROSOLS
14.3. Transport hazard class(es): 2
14.4. Packing group: -
 Hazard label: 2.1



Classification code: 5F
 Limited quantity: 1 L / 30 kg
 Excepted quantity: E0
 Transport category: 2
 Tunnel restriction code: D

Inland waterways transport (ADN)

14.1. UN number or ID number: UN 1950
14.2. UN proper shipping name: AEROSOLS
14.3. Transport hazard class(es): 2
14.4. Packing group: -
 Hazard label: 2.1



Classification code: 5F
 Limited quantity: 1 L / 30 kg
 Excepted quantity: E0

Marine transport (IMDG)

14.1. UN number or ID number: UN 1950
14.2. UN proper shipping name: AEROSOLS (Pentane)
14.3. Transport hazard class(es): 2.1
14.4. Packing group: -
 Hazard label: 2.1



Marine pollutant: Yes
 Limited quantity: 1000 mL / 30 kg
 Excepted quantity: E0
 EmS: F-D, S-U

Air transport (ICAO-TI/IATA-DGR)

14.1. UN number or ID number: UN 1950
14.2. UN proper shipping name: AEROSOLS, FLAMMABLE

Safety Data Sheet

according to UK REACH Regulation

GB600 STERILIT Power Systems Oil spray

Revision date: 20.02.2024

Product code: 00250-0019

Page 9 of 10

14.3. Transport hazard class(es): 2.1**14.4. Packing group:** -

Hazard label: 2.1



Limited quantity Passenger: 30 kg G

Passenger LQ: Y203

Excepted quantity: E0

IATA-packing instructions - Passenger: 203

IATA-max. quantity - Passenger: 75 kg

IATA-packing instructions - Cargo: 203

IATA-max. quantity - Cargo: 150 kg

14.5. Environmental hazards

ENVIRONMENTALLY HAZARDOUS: No

14.6. Special precautions for user

Handle in accordance with good industrial hygiene and safety practice.

14.7. Maritime transport in bulk according to IMO instruments

The transport takes place only in approved and appropriate packaging.

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****EU regulatory information**

Restrictions on use (REACH, annex XVII):

Entry 3, Entry 28, Entry 40

Directive 2004/42/EC on VOC in
paints and varnishes: 95 %Information according to Directive
2012/18/EU (SEVESO III): P3a FLAMMABLE AEROSOLS**National regulatory information**

Employment restrictions: Observe restrictions to employment for juveniles according to the 'juvenile work protection guideline' (94/33/EC). Observe employment restrictions under the Maternity Protection Directive (92/85/EEC) for expectant or nursing mothers.

Water hazard class (D): 2 - obviously hazardous to water

15.2. Chemical safety assessment

For this substance a chemical safety assessment has not been carried out.

SECTION 16: Other information**Changes**

This data sheet contains changes from the previous version in section(s): 9.

Safety Data Sheet

according to UK REACH Regulation

GB600 STERILIT Power Systems Oil spray

Revision date: 20.02.2024

Product code: 00250-0019

Page 10 of 10

Abbreviations and acronyms

ADR = Accord européen relatif au transport international des marchandises Dangereuses par Route
 RID = Règlement concernant le transport international ferroviaire de marchandises dangereuses
 ADN = Accord européen relatif au transport international des marchandises dangereuses par voie de navigation intérieure
 IMDG = International Maritime Code for Dangerous Goods
 IATA/ICAO = International Air Transport Association / International Civil Aviation Organization
 MARPOL = International Convention for the Prevention of Pollution from Ships
 IBC-Code = International Code for the Construction and Equipment of Ships Carrying Dangerous Chemicals in Bulk

GHS = Globally Harmonized System of Classification and Labelling of Chemicals
 REACH = Registration, Evaluation, Authorization and Restriction of Chemicals
 CAS = Chemical Abstract Service
 EN = European norm
 ISO = International Organization for Standardization
 DIN = Deutsche Industrie Norm
 PBT = Persistent Bioaccumulative and Toxic
 vPvB = Very Persistent and very Bio-accumulative

LD = Lethal dose
 LC = Lethal concentration
 EC = Effect concentration
 IC = Median immobilisation concentration or median inhibitory concentration

Classification for mixtures and used evaluation method according to GB CLP Regulation

| Classification | Classification procedure |
|-------------------------|--------------------------|
| Aerosol 1; H222-H229 | On basis of test data |
| Aquatic Chronic 3; H412 | Calculation method |

Relevant H and EUH statements (number and full text)

H220 Extremely flammable gas.
 H222 Extremely flammable aerosol.
 H224 Extremely flammable liquid and vapour.
 H229 Pressurised container: May burst if heated.
 H280 Contains gas under pressure; may explode if heated.
 H304 May be fatal if swallowed and enters airways.
 H336 May cause drowsiness or dizziness.
 H411 Toxic to aquatic life with long lasting effects.
 H412 Harmful to aquatic life with long lasting effects.
 EUH066 Repeated exposure may cause skin dryness or cracking.

Further Information

Data of items 4 to 8, as well as 10 to 12, do partly not refer to the use and the regular employing of the product (in this sense consult information on use and on product), but to liberation of major amounts in case of accidents and irregularities.

The information describes exclusively the safety requirements for the product(s) and is based on the present level of our knowledge.

The delivery specifications are contained in the corresponding product sheet.

This data does not constitute a guarantee for the characteristics of the product(s) as defined by the legal warranty regulations.

(n.a. = not applicable; n.d. = not determined)

(The data for the relevant ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)