



Ultralight & powerful



Air motors



Air motors	AM-25 LT	AM-25 A	AM-25 A BC	AM-25 A RM	AM-25 E
Hose coupling Hose side connection:	Roto Quick with light connection	Roto Quick without light connection	Fixed connection Borden 2(3) hole	Fixed connection standard 4 hole	Roto Quick without light connection
Light:	yes	–	–	–	–
Spray:	internal	internal	internal	internal	external
Speed range [rpm]:	at 2,2 bar: 5,000 – 20,000* at 3 bar: 5,000 – 25,000*				
Rotary direction:	forward/reverse drive				



Air motors	AM-25 E BC	AM-25 E RM	AM-25	AM-25 BC	AM-25 RM
Hose coupling Hose side connection:	Fixed connection Borden 2(3) hole	Fixed connection standard 4 hole	Roto Quick without light connection	Fixed connection Borden 2(3) hole	Fixed connection standard 4 hole
Light:	–	–	–	–	–
Spray:	external	external	–	–	–
Speed range [rpm]:	at 2,2 bar: 5,000 – 20,000* at 3 bar: 5,000 – 25,000*				
Rotary direction:	forward/reverse drive				

* Performance and speed data of air driven motors are largely dependent on the quality of the turbine hoses used and may therefore differ from the specified values.

Order the air motors from your local retail dental supplier.

W&H Dentalwerk Bürmoos GmbH
 Ignaz-Glaser-Straße 53, Postfach 1
 5111 Bürmoos, Austria
 t + 43 6274 6236-0
 f + 43 6274 6236-55
 office@wh.com
 wh.com