



## Thermafil®

Starter Kit

---



Artikelnr.: 216391  
Hersteller: Dentsply Sirona Europe  
Gewicht: 1.28 kg  
Höhe: 110 mm  
Breite: 210 mm  
Länge: 330 mm

**€ 559,00**

*Zzgl. MwSt.*

---

**Starter Kit 1** Thermaprep Ofen 230 V, 1 Thermafil Kit, 20 Obturatoren (ISO 020, ISO 025, ISO 030, ISO 035, ISO 040) Verifier-Kit (6 x Verifier ISO 020-045, 25 mm)