



## onetime perfect

Standardpackung 50 ml Base, 50 ml Katalysator Wash, Zubehör



Artikelnr.: 30223  
Hersteller: DETAX GmbH  
Gewicht: 0.254 kg  
Höhe: 126 mm  
Breite: 58 mm  
Länge: 105 mm

**€ 50,56**

*Zzgl. MwSt.*

**Standardpackung** 50 ml Base, 50 ml Katalysator Wash, 12 Mischkanülen rosa, 12 Intra-Oral-Tips weiß

2-phases, 1-time system on A-silicon basis, for the onetime-technique. Extremely short oral retention (2 min.) due to thermally controlled setting characteristics. Simultaneous application of Putty or Heavy Material and Wash in one working step.

- Putty - very high consistency, kneadable, fast setting
- Heavy - high consistency, fast setting
- Wash - low consistency, free-flowing. Superhydrophilic, thixotropic, extremely sharp definition
- Contrast colors: Putty and Heavy - turquoise, Wash - violet